





Embolectomy Catheters - Single or Double Lumen




ISOMed Embolectomy Catheters - single or double lumen are made with X.R.O. PEBAX®. They are colour coded for easy identification. The latex balloon, fixed with two sutures, is perfectly symmetrical.

Embolectomy catheters are graduated every 10 cm from the distal end, and re-inforced with a removable steel wire. The balloon is protected from sunlight by an opaque sheath.



Embolectomy Catheters - Single or Double Lumen



Embolectomy Catheters Single Lumen

CODE	COLOUR	EXTERNAL Ø	LENGTH	MAX. INFLATE	MAX. BALLOON Ø
300240	Grey	2F (0,70 mm)	40 cm	0,05 ml	4 mm
300260	Grey	2F (0,70 mm)	60 cm	0,05 ml	4 mm
300340	Green	3F (1,00 mm)	40 cm	0,20 ml	5 mm
300380	Green	3F (1,00 mm)	80 cm	0,20 ml	5 mm
300440	Red	4F (1,35 mm)	40 cm	0,75 ml	8 mm
300480	Red	4F (1,35 mm)	80 cm	0,75 ml	8 mm
300540	White	5F (1,70 mm)	40 cm	1,50 ml	10 mm
300580	White	5F (1,70 mm)	80 cm	1,50 ml	10 mm
300640	Blue	6F (2,00 mm)	40 cm	2,00 ml	13 mm
300680	Blue	6F (2,00 mm)	80 cm	2,00 ml	13 mm
300780	Yellow	7F (2,35 mm)	80 cm	2,50 ml	14 mm
300880	Brown	8F (2,70 mm)	80 cm	3,00 ml	15 mm

Other length available upon request : 20, 40, 60, 80, 100 cm

Embolectomy Catheters Double Lumen

CODE	COLOUR	EXTERNAL Ø	LENGTH	MAX. INFLATE	MAX. BALLOON Ø
500480	Red	4 F (1,35 mm)	80 cm	0,75 ml	8 mm
500580	White	5 F (1,70 mm)	80 cm	1,50 ml	10 mm
500680	Blue	6 F (2,00 mm)	80 cm	2,00 ml	13 mm
500780	Yellow	7 F (2,33 mm)	80 cm	2,50 ml	14 mm

Other length available upon request : 20, 40, 60, 80, 100 cm