



Disposable Vein Strippers



CODE	NUMBER OF WIRES	LENGTH	OLIVES 6, 9, 12 AND 15(MM)	HANDLE	PACKAGING (UNITS/BOX)
2004 101	1	100	4	1	10
2004 102	2	100	4	1	
2000 101	1	100	0	0	
2000 102	2	100	0	0	

Sterile double-packaging - Ethylene oxide - Single use



The total flexibility of the nylon wire facilitates the passage across venous valves and minimizes the number of incisions. In addition, it is torsion and kinking free. The smooth wire guarantees easy and traumatic free passage within the veins as well as quick and problem-free shedding of the stripper veins. Thus the extended varicose veins can be quickly controlled for tributary venous vessels and these can be searched specifically and then ligated. This insures a traumatic free and time saving application.

The strippers are available with or without olives and handle,

to allow cost reductions when olives and handle are not used.

A broad handle ensures effective and time-saving stripping of the veins. The handle may be set at either end.

The special designed olives, which can be clipped on each wire end after the wire passes through the vein, enable a retrograde as well as an anterograde stripping, independent from distal or proximal access.

The vein strippers are provided with four different olives of 6, 9, 12 and 15 mm diameter.