A photograph of several petri dishes containing microbial cultures. The top dish shows a greenish-yellow surface, while the bottom dish shows a vibrant red surface. The dishes are arranged in a stack, with the bottom one being the most prominent.

**MICROBIOLOGY
& LIFE SCIENCE
REAGENTS**

PRODUCT LIST

Inspired by knowledge

Index.

Corporate profile	4
--------------------------	---

Microbiology

Culture media.	Dehydrated culture media	9
	Ready-to-use media	18
	Plates 90 mm (petri dish)	18
	Tubes	19
	Flask	20
	Bags	22
	Water plates 55 mm	23
	Contact plates 60 mm (Rodac)	23
	Lock&Block Pharma plates	24
	Hemocultures	24
	Mycobacteria testing	24
	CondaChrome®	25

Supplements and additives.	26
-----------------------------------	----

Bacteriological stains.	27
--------------------------------	----

Antimicrobial susceptibility testing.	Antibiotic discs	28
--	------------------	----

Ingredients.	Agars	30
	Peptones	30
	Carbohydrates	30

Molecular biology

Agaroses.	Standard agaroses	34
	Agarose for DNA recovery	34
	Agarose for high DNA resolution	34
	Special agaroses	34
	Agarose reagents	34
	Low melting agaroses	35
	Molecular screening	35

Selection antibiotics.	35
-------------------------------	----

Microbiological culture media.	36
---------------------------------------	----

PCR - Real Time detection kits

Condagene®	Extraction kits	39
	Pathogen detection kits	39

Cell culture

Condalow®	Condalow®	41
------------------	-----------	----


Additional information.	42
--------------------------------	----

Who we are

A leading European manufacturer of culture media

Founded in 1960 we are one of the main European manufacturers of dehydrated culture media for microbiology and molecular biology, which positions us today as a leading international privately held company.

From our factory in Madrid, Spain we export to more than 130 countries all over the world, directly and through authorized distributors. Our wide distribution network along with our highly skilled staff, quality and specialized product range, are the keys to our success.



Experts.
Mastering culture media.

Flexible & reliable.
Tailored to your needs.

Innovative.
Inspiring the future.


What we do

Develop, manufacture and distribute high quality culture media for micro and molecular biology

The experience acquired in the design and manufacture of culture media has matured us as specialists in the field of microbiology. We design and produce more than 700 dehydrated culture media. Condalab is also recognized for providing key ingredients such as agars, peptones and agaroses, among other products. Our catalog also includes media solutions for a wide range of molecular biology application



Microbiology.
Dehydrated culture media
Ready to use culture media
Supplements
Antimicrobial susceptibility testing
Stains



Molecular biology.
Dehydrated culture media
Agaroses
Molecular biology stains



Bioingredients
Agars
Peptones
Carbohydrates

Who we do it for

We have forged an increasingly large network of clients as a result of our ability to listen as well as our commitment towards finding optimal solutions.

Condalab products address the following sectors:



Clinical analyses.

Hospitals, veterinary clinics, clinical and food analysis laboratories.



Quality controls.

Food, beverage, breweries. Pharmaceutical industry, veterinary and cosmetic products.



Production processes.

Fermentation processes, vaccine manufacturing, probiotic products, prepared media manufacturers.



R&D.

Laboratories, research centers, universities.



We follow strict quality control procedures before, during and after the manufacturing process.

How we do it

We take compliance seriously

Condalab continuously improves and increases its production to meet the highest international quality standards. We comply with the ISO 9001:2015, ISO 13485:2018 and European Community CE marking for in vitro medical devices.

Our formulations comply with various international standards such as ISO, European Pharmacopoeia, FDA, APHA, USP and AOAC. We follow strict quality control procedures before, during and after the manufacturing process to guarantee products of absolute quality and consistency from batch to batch.

MICROBIOLOGY



Culture media.

Endorsed by more than 60 years of experience in the design and manufacturing of products for micro and molecular biology, Condalab offers one of the widest ranges of microbiological culture media on the market. The high quality raw materials used in its production, as well as the strict quality controls carried out throughout the process, makes Condalab's culture media the best option for your laboratory; guaranteeing maximum recovery and performance. The product catalog includes more than 700 media in dehydrated format as well as ready-to-use.

Dehydrated culture media.

Dehydrated culture media has a wide variety of applications in industrial microbiology (food industry, cosmetics, pharmaceuticals, vaccine production, animal health, analysis laboratories, etc.), as well as clinical microbiology (human and veterinary). The dehydrated powder form gives the culture medium a longer lifetime of up to 5 years, allowing them to be prepared easily.

CAT.	PRODUCT	PACK SIZE
1252	A1 Medium	500 grams
1391	Acetamide Agar	500 grams
1211	Acetamide Broth	500 grams
1155	Acetamide Broth Base ISO	500 grams
2017	Acetamide Broth ISO	500 grams
1192	Acetate Differential Agar	500 grams
1459	Actinomycete Isolation Agar with Glycerol	500 grams
1370	Agar Base Aeromonas (RYAN)	500 grams
2022	Agar Test pH 6,0 for Inhibition Test	500 grams
2024	Agar Test pH 8,0 for Inhibition Test	500 grams
2002	Alicyclobacillus Detection Agar	500 grams
1407	Alkaline Peptone Water	500 grams
2155	Alkaline Peptone Water ISO	500 grams
1530	Amies Transport Medium without Charcoal	500 grams
1535	Amies Transport with Charcoal	500 grams
1000	Anaerobic Agar	500 grams
1520	Antibiotic Medium N° 1 (Seed Agar) USP	500 grams
1528	Antibiotic Medium N° 11 (Neomycin Test)	500 grams
1525	Antibiotic Medium N° 12	500 grams
1386	Antibiotic Medium N° 23	500 grams
1524	Antibiotic Medium N° 5 (Streptomycin test) USP	500 grams
1004	Antibiotic Medium N° 8 USP	500 grams
1537	Antibiotic Medium N°10 EP	500 grams
1002	Antibiotic Medium N°2 (Base Agar)	500 grams
1534	Antibiotic Medium N°3 USP	500 grams
2049	APT Agar (All Purpose Tween)	500 grams
2067	APT Broth	500 grams
1207	Asparagine Broth	500 grams
1113	Azide Blood Agar Base	500 grams
1422	Azide Dextrose Broth	500 grams
1124	Bacillus Cereus Selective Agar Base (MYP)	500 grams
1343	Bacillus Cereus Selective Agar Base (MYP) ISO	500 grams
1100	Baird Parker Agar Base	500 grams
2053	BAT Broth	500 grams
1051	BCP Agar	500 grams
1320	BCP Glucose Agar	500 grams

CAT.	PRODUCT	PACK SIZE
2009	BCP Glucose Broth	500 grams
1006	Biggy Agar	500 grams
1031	Bile Esculin Agar ISO	500 grams
1372	Bile Esculin Azide Agar	500 grams
1005	Bile Esculin Azide Agar ISO	500 grams
1359	Bile Esculin Azide Broth	500 grams
1011	Bismuth Sulfite Agar (Wilson Blair) USP	500 grams
1108	Blood Agar Base	500 grams
1128	Blood Agar Base + Nalidixic Acid	500 grams
1119	Blood Agar Base Low pH	500 grams
1328	Blood Agar Base N°2 ISO	500 grams
1441	Bolton Selective Enrichment Broth Base ISO	500 grams
1107	Bordet-Gengou Agar Base	500 grams
1490	Bordet-Gengou Agar with Charcoal	500 grams
1451	BPRM Broth Base (Bacteroides Phage Recovery Medium)	500 grams
1048	Brain Heart Infusion Agar (BHI Agar)	500 grams
1400	Brain Heart Infusion Broth (BHI Broth)	500 grams
1331	Brain Heart Infusion Broth ISO	500 grams
2006	Brettanomyces Agar	500 grams
2005	Brettanomyces Selective Broth	500 grams
1078	Brilliant Green Agar	500 grams
1143	Brilliant Green Agar Modified ISO	500 grams
1366	Brilliant Green Agar with Sulfadiazine	500 grams
1010	Brilliant Green Bile Agar	500 grams
1228	Brilliant Green Bile Broth 2% ISO	500 grams
1253	Brilliant Green Tetrathionate Bile Broth	500 grams
1012	Brucella Agar	500 grams
1223	Brucella Broth	500 grams
1374	Brucella Medium Base	500 grams
1247	Bryant-Burkey Broth Base (modif. with resazurin)	500 grams
1277	Buffered Listeria Enrichment Broth with Supplements	500 grams
1401	Buffered Peptone Water EP/USP	500 grams
1402	Buffered Peptone Water ISO	500 grams
1406	Buffered Saline Peptone Water	500 grams
1347	Burkholderia Cepacia Agar Base	500 grams
2201	Burkholderia cepacia Selective Agar Base USP	500 grams
1184	Buti Medium	500 grams
1069	Calcium Caseinate Agar	500 grams
1131	Campylobacter Agar Base (Preston)	500 grams
1129	Campylobacter Agar Base Blood Free (CCDA) ISO	500 grams
2166	Campylobacter Preston Broth Base ISO	500 grams
1342	Carbohydrates Utilization Broth Base ISO	500 grams
1529	Cary-Blair Medium	500 grams
1102	Cetrimide Agar Base EP/USP/ISO	500 grams
1017	Chapman Stone Agar	500 grams
1301	Chloramphenicol Agar (YGC Agar) ISO	500 grams
1016	CLED Agar (Cystine Lactose Electrolyte Deficient)	500 grams
1303	CLED Agar with Andrade's Indicator	500 grams
1447	Clostridium Difficile Agar Base	500 grams
1132	Clostridium Perfringens Agar Base (m-CP)	500 grams

CAT.	PRODUCT	PACK SIZE
1104	Columbia Agar Base EP/USP/ISO	500 grams
1229	Columbia Broth	500 grams
1152	Columbia CNA Agar Base	500 grams
2222	Cooked Meat Medium	500 grams
1164	Corn Meal Agar	500 grams
2143	Cronobacter Selective Broth (CSB) ISO	500 grams
1502	CTA Medium	500 grams
1250	Czapek-Dok Modified Broth	500 grams
1015	Czapek-Dox Modified Agar	500 grams
1045	DCLS Agar (Desoxycholate, Citrate, Lactose, Sucrose)	500 grams
1067	Deoxycholate Citrate Agar	500 grams
1157	Dermatophytes Agar Base	500 grams
1020	Desoxycholate Agar	500 grams
1025	Desoxycholate Lactose Agar	500 grams
1421	DEV Lactose Broth	500 grams
1021	Dextrose Agar	500 grams
1203	Dextrose Broth (Glucose Broth)	500 grams
1146	Dextrose Tryptone Agar with Starch	500 grams
2003	Dey-Engley Neutralizing Broth ISO	500 grams
1185	DHL Agar	500 grams
1161	Dichloran Glycerol Agar (DG 18) ISO	500 grams
1416	Differential Reinforced Clostridium Broth	500 grams
1028	DNase Test Agar (Deoxyribonuclease Activity)	500 grams
2028	DNase Test Agar with Toluidine Blue (Deoxyribonuclease Activity)	500 grams
1337	E. coli O157:H7 Presumptive ID Agar	500 grams
1522	EC Medium ISO	500 grams
1180	EC with MUG Fluorogenic Agar	500 grams
1362	EE Broth	500 grams
1539	Elliker Medium	500 grams
1118	Endo Agar Base	500 grams
1276	Endo Base Broth	500 grams
1137	Endo LES Agar Base	500 grams
1018	Enterococcus Confirmatory Agar	500 grams
1070	Enterococcus Selective Agar (Enterococcosel Agar)	500 grams
1204	Enterococcus Selective Broth (Enterococcosel Broth)	500 grams
1039	Eosin Methylene Blue Agar (EMB)	500 grams
1574	Esculin Iron Agar	500 grams
1254	ESTY Broth	500 grams
1555	ESTY Medium	500 grams
1036	Eugon Agar	500 grams
2152	Eugon LT 100 Agar	500 grams
2110	Eugon LT 100 Broth	500 grams
1230	EVA Broth (Ethyl, Violet, Azide, Litsky)	500 grams
1212	Ewing Malonate Broth Modified	500 grams
1127	Fecal Coliforms Agar Base (m-FC)	500 grams
1121	Fecal Coliforms Broth Base (m-FC)	500 grams
1264	FMM Broth	500 grams
1106	GC Agar Base	500 grams
2075	GC Broth	500 grams
1526	Gelatin Lactose Medium	500 grams

CAT.	PRODUCT	PACK SIZE
1187	Gelysate Agar (Plate Count Agar Sugar Free)	500 grams
1232	Giolitti-Cantoni Broth	500 grams
1287	Giolitti-Cantoni Broth ISO	500 grams
1094	Glucose Chloramphenicol Agar	500 grams
1258	Glucose Chloramphenicol Broth	500 grams
2150	Glucose OF Medium ISO	500 grams
1248	GN Enrichment Broth (Hajna)	500 grams
2036	Granada Agar Base	500 grams
2026	GYC Agar	500 grams
1464	Haemophilus Test Medium	500 grams
2001	HC Agar Base	500 grams
1323	Heart Infusion Broth	500 grams
1030	Hektoen Enteric Agar ISO	500 grams
2066	Helicobacter Agar Modified	500 grams
1165	Hoyle Medium Base	500 grams
1504	Indole Nitrite Medium (Trypticasein Nitrate Medium)	500 grams
2097	Inositol Brilliant Green Bile Agar (Plesiomonas Differential Agar)	500 grams
1361	Irgasan Ticarcillin and Potassium Chlorate Broth (ITC) ISO	500 grams
1085	Iron Agar (Lyngby)	500 grams
1559	Iron Sulfite Agar ISO	500 grams
1001	Iso-Sensitest Agar	500 grams
1027	KAA Confirmatory Agar	500 grams
1209	KAA Presumptive Broth	500 grams
1034	KF Streptococcal Agar	500 grams
1101	KF Streptococcal Agar with Bromocresol Purple	500 grams
1346	KF Streptococcal Broth	500 grams
1531	King A Medium (Pseudomonas P) Agar USP/ISO	500 grams
1532	King B Medium (Pseudomonas F Agar) ISO	500 grams
1053	King FG Agar	500 grams
1042	Kliger Iron Agar	500 grams
1364	Kligler Iron Agar	500 grams
1200	Koser Citrate Broth	500 grams
1344	Lactose Agar with Bromothymol Blue and Crystal Violet (Drigalski)	500 grams
1206	Lactose Broth	500 grams
1272	Lactose Broth with Bromocresol Purple	500 grams
2182	Lactose Glucuronide m-Agar (MLGA)	500 grams
1420	Lactose Peptone Broth DEV	500 grams
1009	Lactose Sulfite Broth Base ISO	500 grams
1309	Lauryl Sulfate Agar	500 grams
1310	Lauryl Sulfate Broth (Lauryl Tryptose Broth – LTB) ISO	500 grams
1445	Lauryl Sulphate Tryptose Broth Modified (mLST)	500 grams
1311	Legionella BCYE Agar Base ISO	500 grams
1111	Lethen Agar Modified	500 grams
1244	Lethen Broth Modified ISO/BAM	500 grams
1050	Levine Agar (EMB) ISO/BAM	500 grams
1389	LIM Medium (Lysine, Indole, Motility Medium)	500 grams
2068	Linden Grain Medium	500 grams
1133	Listeria Agar Base Oxford ISO	500 grams
1141	Listeria Agar Base Palcam ISO	500 grams
1120	Listeria Enrichment Broth Base ISO	500 grams

CAT.	PRODUCT	PACK SIZE
1077	Listeria Enrichment Broth With Supplements	500 grams
1182	Listeria Fraser Broth Base ISO	500 grams
1183	Listeria Half-Fraser Broth Base ISO	500 grams
1142	Liver Agar	500 grams
1242	Liver Broth	500 grams
2208	Liver Veal Agar	500 grams
1084	LMDA Agar	500 grams
2087	Loeffler Serum Medium Base	500 grams
2149	L-Ornithine Decarboxylation Medium ISO	500 grams
1439	Lowenstein Jensen Medium Base w/o M.Green	500 grams
1399	LPT Dilution Broth	500 grams
1208	Lysine Decarboxylase Broth	500 grams
1176	Lysine Decarboxylase Medium ISO	500 grams
1044	Lysine Iron Agar	500 grams
1318	M17 Agar	500 grams
1052	Macconkey Agar EP/USP/ISO	500 grams
1035	MacConkey Agar N° 2	500 grams
1098	MacConkey Agar w/o Crystal Violer and w/o Sodium Chloride	500 grams
1099	Macconkey Agar with Sorbitol (CT-SMAC) ISO	500 grams
1210	MacConkey Broth EP/USP	500 grams
1573	MAL Agar	500 grams
1275	Malonate Phenylalanine Broth	500 grams
1038	Malt Extract Agar	500 grams
1245	Malt Extract Broth	500 grams
1509	Mannitol Nitrate Motility Medium	500 grams
1062	Mannitol Salt Agar (MSA) (Chapman Medium) EP/USP/ISO	500 grams
1487	Mannitol Salt Broth	500 grams
1059	Marine Agar	500 grams
1217	Marine Broth	500 grams
2133	m-Dixon agar	500 grams
2108	M-Green Yeast and Mold Agar	500 grams
1365	Minerals Modified Glutamate Broth (MMBG) ISO	500 grams
1510	MIO Medium (Motility-Indole-Ornithine)	500 grams
1393	MLO Medium	500 grams
1438	Modified Nocive Brewers Bacteria Agar Base	500 grams
1186	Modified Salmonella Shigella Agar	500 grams
1452	Modified Scholten's Broth	500 grams
1376	Modified Semisolid Rappaport Vassiliadis Medium (MSRV) ISO	500 grams
1112	Moeller KCN Broth Base	500 grams
1202	Mossel EE Broth EP/USP	500 grams
2116	Motility Test Medium	500 grams
1043	MRS Agar	500 grams
1433	MRS Agar Low pH ISO	500 grams
1215	MRS Broth	500 grams
1431	MRS Broth Low pH	500 grams
1295	MRS Broth w/o Dextrose	500 grams
1288	MRS Broth w/o Dextrose and w/o Beef Extract	500 grams
1512	MR-VP Medium	500 grams
2216	m-TEC Agar	500 grams
1058	Mueller Hinton Agar	500 grams

CAT.	PRODUCT	PACK SIZE
1055	Mueller Hinton Agar II	500 grams
1214	Mueller Hinton Broth	500 grams
1130	Muller Kauffmann Tetrathionate Broth Base	500 grams
1173	Muller Kaufmann Broth Base with brilliant Green and Novobiocin (MKTTN) ISO	500 grams
1072	Mycobiotic Agar (Fungal Selective Agar)	500 grams
1449	Neutralizing Agar	500 grams
1281	Nitrate Broth	500 grams
1565	Nitrate Motility Base Medium	500 grams
1440	Nocive Brewers Bacteria Broth Base Modified	500 grams
1314	Nutrient Agar Dev Regulations	500 grams
1060	Nutrient Agar ISO	500 grams
1156	Nutrient Agar ISO	500 grams
1355	Nutrient Agar with Sodium Chloride ISO	500 grams
1216	Nutrient Broth	500 grams
2107	Nutrient Broth	500 grams
1251	Nutrient Broth N°2	500 grams
1300	Nutrient Gelatin	500 grams
2060	Oatmeal Agar	500 grams
1500	OF Basal Medium (Hugh And Leifson)	500 grams
1527	OGA Medium (Oxytetracycline Glucose Agar Base) (OGYE)	500 grams
1307	Orange Serum Agar	500 grams
1057	Osmophilic Agar	500 grams
1403	Peptone Water (Tryptone Water) ISO	500 grams
1357	Peptone Water with Lactose	500 grams
1115	Phenol Red Broth Base	500 grams
1023	Phenol Red Dextrose Agar	500 grams
1235	Phenol Red Dextrose Broth	500 grams
1239	Phenol Red Sucrose Broth	500 grams
1040	Phenylalanine Agar	500 grams
1022	Potato Dextrose Agar EP/USP	500 grams
1261	Potato Dextrose Broth	500 grams
1140	PPLO Agar Base w/o Crystal Violet	500 grams
1262	PPLO Broth Base w/o Crystal Violet	500 grams
2061	Presence-Absence Broth	500 grams
1356	Pseudomonas CFC Agar Base ISO	500 grams
1153	Pseudomonas CN Agar Base ISO	500 grams
1071	R2A Agar EP/USP	500 grams
1470	R2A Agar w/o Agar	500 grams
1061	Raka-Ray Agar Base	500 grams
1174	Rappaport Soy Broth (Vassiliadis) ISO	500 grams
1414	Rappaport Vassiliadis Broth EP/USP	500 grams
1087	Reinforced Clostridal Agar	500 grams
1007	Reinforced Clostridal Medium EP/USP	500 grams
1096	Rogosa SL Agar	500 grams
1234	Rogosa SL Broth	500 grams
1081	Rose Bengal Agar + Chloramphenicol	500 grams
1160	Rose Bengal Agar + Cholamphenicol + Dichloran (DRBC Agar) ISO	500 grams
2055	Rosenow Broth	500 grams
2035	Rothe Agar (Glucose Agar w/ Azide)	500 grams
1238	Rothe Broth (Glucose Broth With Azide)	500 grams

CAT.	PRODUCT	PACK SIZE
1090	Sabouraud Dextrose Agar + Chloramphenicol	500 grams
1089	Sabouraud Dextrose Agar + Chloramphenicol + Cycloheximide	500 grams
2059	Sabouraud Dextrose Agar + Chloramphenicol + Gentamicin	500 grams
1134	Sabouraud Dextrose Agar + Chloramphenicol EP/USP/ISO	500 grams
1088	Sabouraud Dextrose Agar + Cycloheximide	500 grams
1166	Sabouraud Dextrose Agar 2%	500 grams
1024	Sabouraud Dextrose Agar EP/USP/ISO	500 grams
1284	Sabouraud Dextrose Broth + Chloramphenicol	500 grams
1205	Sabouraud Dextrose Broth EP/USP	500 grams
1506	Sabouraud Fluid Medium	500 grams
1379	Saline Lactose Broth with Phenol Red	500 grams
1405	Saline Peptone Water ISO	500 grams
1064	Salmonella Shigella Agar (SS Agar)	500 grams
1360	Salmonella Shigella Agar with Sodium Deoxycholate and Calcium Chloride (SSDC)	500 grams
1463	Salmonella Shigella Broth	500 grams
1326	Salt Egg Yolk Agar Base	500 grams
1413	San Francisco Base Medium Modified	500 grams
2191	SCDLP 80 Broth ISO	500 grams
1066	Schaedler Agar	500 grams
1218	Schaedler Broth	500 grams
1282	Schubert Broth	500 grams
1220	Selenite Cystine Broth	500 grams
1065	Sellers Agar	500 grams
2046	Semisolid Nutrient Agar ISO	500 grams
2078	Shigella broth ISO	500 grams
1514	SIM Medium	500 grams
1014	Simmons Citrate Agar ISO	500 grams
1109	Slanetz-Bartley Medium ISO	500 grams
1435	Slanetz-Bartley Medium w/o TTC	500 grams
1158	Sodium Chloride Buffered Peptone Solution EP	500 grams
1222	Sodium Selenite Broth	500 grams
1298	Sorbitol Peptone Broth And Bile Salts (PSB) ISO	500 grams
1082	SPS Agar (Sulfite Polymyxin Sulfadiazine)	500 grams
1178	Standard Count Agar	500 grams
1056	Standard Methods Agar (PCA) ISO/APHA	500 grams
1033	Standard Methods Agar with Powdered Milk APHA/ISO	500 grams
1177	Standard Nutrient Agar I	500 grams
1286	Standard Nutrient Broth I	500 grams
1032	Staphylococcus Agar N° 110	500 grams
2099	Streptococcus Selective Agar	500 grams
2120	Streptococcus Selective Agar Modified (without Neomycin without Colistin)	500 grams
1518	Stuart Transport Medium	500 grams
1315	Syncase Broth	500 grams
1029	T.S.C. Agar Base (Tryptose Sulfite Cycloserine) ISO	500 grams
1397	TAT Broth Base (CPLP)	500 grams
1074	TCBS Agar ISO	500 grams
1114	Tetrathionate Broth Base 500 g	500 grams
1190	TGEA Medium (Tryptone Glucose Yeast Extract Agar)	500 grams
1241	Thioglycollate Broth NIH	500 grams
1508	Thioglycollate Fluid Medium EP/USP	500 grams

CAT.	PRODUCT	PACK SIZE
1516	Thioglycollate Fluid Medium without Indicator	500 grams
1533	Thioglycollate Medium ISO	500 grams
2029	Thrichomonas Broth (Kupferberg)	500 grams
1317	Tinsdale Agar Base	500 grams
1236	Todd Hewitt Broth	500 grams
2011	TOS Propionate Agar Base ISO	500 grams
1591	TPY Broth	500 grams
1046	Triple Sugar Iron Agar (TSI)	500 grams
1172	Triple Sugar Iron Agar (TSI) ISO	500 grams
1041	Trypticasein Glucose Extract Agar	500 grams
1068	Trypticasein Soy Agar (TSA) EP/USP/ISO	500 grams
1561	Trypticasein Soy Agar (TSA) N°2	500 grams
1189	Trypticasein Soy Blood Agar Base	500 grams
1224	Trypticasein Soy Broth (TSB) EP/USP/ISO	500 grams
2157	Trypticasein Soy Broth (TSB) With Tween 80	500 grams
1292	Trypticasein Soy Broth Modified with Novobiocin (mTSB) ISO	500 grams
1225	Trypticasein Soy Broth w/o Dextrose	500 grams
1013	Tryptone Bile Salts Agar	500 grams
1237	Tryptophan Culture Broth ISO	500 grams
1047	Tryptose Agar	500 grams
1162	Tryptose Blood Agar Base	500 grams
1322	Tryptose Broth	500 grams
1243	Tryptose Phosphate Broth	500 grams
1294	Tryptose Phosphate Broth II	500 grams
1138	Tryptose Soy Agar ISO	500 grams
1581	Tryptose Agar w/o Thiamine	500 grams
1198	TSA Agar N°2 Modified	500 grams
1396	TSA Lethen Agar	500 grams
1075	TSN Agar (Tryptone Sulfite Neomycin)	500 grams
1398	TSYEA Agar (Tryptone Soy Yeast Extract Agar) ISO	500 grams
1339	TSYEB Broth (Tryptone Soy Yeast Extract Broth) ISO	500 grams
1076	TTC Chapman Agar (Lactose Agar TTC with Tergitol 7)	500 grams
1562	Universal Beer Agar (UBA)	500 grams
1110	Urea Agar Base (Christensen) ISO	500 grams
2180	Urea Agar Base (Christensen) ISO	500 grams
1226	Urea Broth	500 grams
1227	Urea Indole Broth ISO	500 grams
1279	UVM-I Listeria Selective Enrichment Broth Modified	500 grams
1280	UVM-II Listeria Selective Enrichment Broth Modified	500 grams
1410	Vancomycin Screen Agar	500 grams
1380	Vegetable Peptone Broth (TSB Vegetable)	500 grams
1092	Violet Red Bile Agar With Glucose (VRBG) EP/USP/ISO	500 grams
1144	Violet Red Bile Agar With Glucose and Lactose (VRBLG) USP	500 grams
1093	Violet Red Bile Agar With Lactose (VRBL) ISO	500 grams
1485	VL Anaerobic Blood Agar	500 grams
1079	Vogel-Johnson Agar	500 grams
1392	VP Semisolid Medium	500 grams
1503	Wilkins Chalgren Medium	500 grams
1568	Wilkins Chalgren Medium II	500 grams
1026	WL Differential Agar	500 grams

CAT.	PRODUCT	PACK SIZE
1086	WL Nutrient Agar	500 grams
1572	WL Nutrient Broth	500 grams
1316	Wort Agar	500 grams
1444	Wort Broth	500 grams
1080	XLD Agar (Xylose Lysine Desoxycholate Agar) EP/USP	500 grams
1274	XLD Agar (Xylose Lysine Desoxycholate Agar) ISO	500 grams
1159	XLT4 Agar Base	500 grams
1049	Yeast Extract Agar (YEA) ISO	500 grams
1312	Yeast Extract Agar for Molds	500 grams
1097	Yeast Extract Soy Agar	500 grams
1268	Yeast Extract Soy Broth (TSB no animal)	500 grams
1194	Yeast Mold Agar	500 grams
2008	Yeast Mold Broth	500 grams
1126	Yersinia Selective Agar Base (CIN) ISO	500 grams
2100	YSG Broth	500 grams

Ready-to-use media.

Culture media preparation requires expert staff and time, in addition to the quality controls of the batch itself. Ready to use culture medium offers industry professionals security, easy use and time optimization, simplifying day-to-day work. Condalab products are available in a wide variety of presentations (90 mm, 55 mm and contact plates) and different medium formats such as tubes, bottles and bags, as well as products for the control of clean rooms and sterile environments.

Plates 90 mm (Petri dish)

CAT.	PRODUCT	PACK SIZE
945	Bacillus Cereus Selective Agar (MYP)	Pack of 20 plates
921	Baird-Parker Agar ISO	Pack of 20 plates
984	Baird-Parker RPF Agar ISO	Pack of 20 plates
920	Biggy Agar	Pack of 20 plates
952	Bile Esculin Azide Agar	Pack of 20 plates
912	Blood Agar Base N° 2 ISO	Pack of 20 plates
937	Brain Heart Infusion Agar (BHI Agar)	Pack of 20 plates
915	Brilliant Green Agar (BGA) ISO	Pack of 20 plates
983	Burkholderia Cepacia Agar (BCSA)	Pack of 20 plates
981	CCA Coliforms Chromogenic Agar	Pack of 20 plates
878	CCDA Campylobacter Agar ISO	Pack of 20 plates
916	Cetrimide Agar EP/USP/ISO	Pack of 20 plates
885	Chocolate Agar with VX	Pack of 20 plates
866	Chromogenic Cronobacter Isolation Agar (CCI) (ISO)	Pack of 20 plates
905	CLED Agar (Cystine Lactose Electrolyte Deficient)	Pack of 20 plates
974	Columbia CNA Agar with 5% Sheep Blood	Pack of 20 plates
931	Columbia with 5% Sheep Blood Agar	Pack of 20 plates
879	Dermatophyte Agar	Pack of 20 plates
907	Dextrose Sabouraud Agar EP/USP	Pack of 20 plates
832	DG 18 Metal Agar	Pack of 20 plates
833	Dichloran Rose Bengal Chloramphenicol Metal Agar	Pack of 20 plates
953	Eugon LT Modified Agar ISO	Pack of 20 plates
812	Gelysate Agar (Plate Count Agar Sugar Free)	Pack of 20 plates
918	Hektoen Enteric Agar ISO	Pack of 20 plates
985	Legionella Agar (CYE / GVPC) ISO	Pack of 20 plates
828	Legionella BCYE with AB Agar ISO	Pack of 20 plates
944	Legionella CYE Agar ISO	Pack of 20 plates
978	Legionella CYE Agar w/o Cysteine ISO	Pack of 20 plates
906	Levine Agar (EMB) ISO	Pack of 20 plates
884	Listeria Chromogenic Agar Base according to Ottaviani and Agosti (ALOA) ISO	Pack of 120 plates
939	Listeria Chromogenic Agar Base according to Ottaviani and Agosti (ALOA) ISO	Pack of 20 plates
996	Luria Agar (Miller's LB Agar)	Pack of 20 plates
997	Luria Agar (Miller's LB Agar) with Ampicillin (100µg/ml)	Pack of 20 plates
998	Luria Agar (Miller's LB Agar) with Kanamycin (50µg/ml)	Pack of 20 plates
900	MacConkey Agar EP/USP/ISO	Pack of 20 plates
827	Malt Extract Agar	Pack of 20 plates
821	Malt Extract Agar Acidified (pH 4,3)	Pack of 20 plates
831	Malt Extract with Chloramphenicol Agar	Pack of 20 plates
917	Mannitol Salt Agar (MSA) (Chapman Medium) EP/USP/ISO	Pack of 20 plates
886	Marine Agar with Cycloheximide	Pack of 20 plates
988	M-Green Yeast and Mould Agar	Pack of 20 plates

CAT.	PRODUCT	PACK SIZE
943	MRS Agar	Pack of 20 plates
813	MRS Agar Low pH ISO	Pack of 20 plates
901	Mueller Hinton Agar	Pack of 20 plates
969	Mueller Hinton Agar with 5% Blood	Pack of 20 plates
902	Nutrient Agar EP/USP/ISO	Pack of 20 plates
934	Nutrient Agar ISO	Pack of 20 plates
851	OGYE Agar ISO	Pack of 20 plates
955	Palcam Listeria Agar ISO	Pack of 20 plates
849	Potato Dextrose Agar EP/USP	Pack of 20 plates
947	Potato Dextrose Agar with Chloramphenicol ISO	Pack of 20 plates
829	Pseudomonas CFC Agar ISO	Pack of 20 plates
993	R2A Agar EP (Irradiated)	Pack of 20 plates
950	Rose Bengal Agar with Chloramphenicol	Pack of 20 plates
843	Sabouraud Chloramphenicol Dextrose Agar. Double Wrapped. Irradiated EP/USP/ISO	Pack of 20 plates
842	Sabouraud Chloramphenicol Dextrose Agar. Triple Wrapped. Irradiated. EP/USP/ISO	Pack of 20 plates
908	Sabouraud Dextrose Agar with Chloramphenicol EP/USP/ISO	Pack of 20 plates
841	Sabouraud Dextrose Agar. Double Wrapped. Irradiated EP/USP/ISO	Pack of 20 plates
838	Sabouraud TLHTh Agar. EP/USP Double Wrapp. Irradiated	Pack of 20 plates
909	Salmonella Shigella Agar (SS Agar)	Pack of 20 plates
946	Slanetz-Bartley Medium ISO	Pack of 20 plates
903	Standard Methods Agar (PCA) ISO/APHA	Pack of 20 plates
982	TBX Chromogenic Agar (Tryptone Bile X-Glucuronide) ISO	Pack of 20 plates
957	TCBS Agar ISO	Pack of 20 plates
904	Trypticasein Soy Agar (TSA) EP/USP/ISO	Pack of 20 plates
836	Trypticasein Soy Agar (TSA) EP/USP/ISO. Doble wrapped. Irradiated	Pack of 20 plates
839	Trypticasein Soy Agar (TSA) EP/USP/ISO. Triple Wrapped. Irradiated	Pack of 20 plates
970	Trypticasein Soy Agar with 5% Sheep Blood	Pack of 20 plates
840	Tryptone Soya Agar (TSA) with Neutralizings. Double Wrapped. Irradiated	Pack of 20 plates
837	TSA TLHTh. (Trypticasein Soy Agar) EP/USP Triple Wrapped Irradiated	Pack of 20 plates
964	TSC Agar with Egg Yolk	Pack of 20 plates
911	Violet Red Bile Agar With Glucose (VRBG) EP/USP/ISO	Pack of 20 plates
910	Violet Red Bile Lactose Agar (VRBL) ISO	Pack of 20 plates
930	XLD Agar (Xylose Lysine Desoxycholate Agar) EP/USP	Pack of 20 plates
847	XLD Agar (Xylose Lysine Desoxycholate Agar) ISO	Pack of 20 plates
941	Yeast Extract Agar (YEA) ISO	Pack of 20 plates
834	Yeast Extract Peptone Dextrose Agar (YPD Agar)	Pack of 20 plates
933	Yersinia CIN Agar ISO	Pack of 20 plates
830	YGC Agar ISO	Pack of 20 plates

Tubes

CAT.	PRODUCT	PACK SIZE
4018	Alkaline Peptone Water ISO	Pack of 20 Tubes
4202	Brain Heart Infusion Broth ISO	Pack of 20 Tubes
4253	Brilliant Green Bile Broth 2% (Doble Concentration) ISO	Pack of 20 Tubes
4213	Brilliant Green Broth 2% ISO	Pack of 20 Tubes
4035	Buffered Peptone Water 0,1% EP/USP/ISO	Pack of 20 Tubes
4037	Buffered Peptone Water 0,1% EP/USP/ISO	Pack of 20 Tubes
4250	Buffered Peptone Water ISO	Pack of 20 Tubes
4127	Coletsos Medium	Pack of 20 Tubes

CAT.	PRODUCT	PACK SIZE
4038	Diluent of Beerens	Pack of 20 Tubes
4265	EC Medium ISO	Pack of 20 Tubes
4046	Eugon LT 100 Broth	Pack of 20 Tubes
4116	Fluid A EP/USP	Pack of 20 Tubes
4261	Giolitti-Cantoni Broth (Double Concentration) ISO	Pack of 20 Tubes
4228	Giolitti-Cantoni Broth ISO	Pack of 20 Tubes
4064	Glucose OF	Pack of 20 Tubes
4026	Glucose OF Medium ISO	Pack of 20 Tubes
4002	Kligler Iron Agar	Pack of 20 Tubes
4208	Lactose Broth	Pack of 20 Tubes
4013	Lactose Broth (Double Concentration)	Pack of 20 Tubes
4135	Lactose Sulphite Broth (C+) ISO	Pack of 20 Tubes
4040	Lauryl Sulfate Broth (Lauryl Tryptose Broth – LTB) C+ ISO	Pack of 20 Tubes
4039	Lauryl Sulfate Broth (Lauryl Tryptose Broth – LTB) ISO	Pack of 20 Tubes
4031	Lethen Broth Modified	Pack of 20 Tubes
4024	Listeria enrichment Broth Fraser ISO	Pack of 20 Tubes
4053	Listeria enrichment Broth Hslf Fraser ISO	Pack of 20 Tubes
4229	Lysine Decarboxylase Broth ISO	Pack of 20 Tubes
4007	Lysine Iron Agar	Pack of 20 Tubes
4047	MacConkey Broth EP/USP	Pack of 20 Tubes
4071	Malt Extract Agar Acidified (pH 4,3)	Pack of 20 Tubes
4103	Malt Extract Broth Acidified (pH 4,3)	Pack of 50 Tubes

Flasks

CAT.	PRODUCT	PACK SIZE
4685	Alkaline Peptone Water ISO	Pack of 10 x 225 ml Flask
4669	Baird Parker RPF Base Agar ISO	Pack of 10 x 90 ml Flask
4677	Baird-Parker Medium Base ISO	Pack of 10 x 90 ml Flask
4678	Baird-Parker Medium Base ISO	Pack of 10 x 180 ml Flask
4658	Bolton Broth Base ISO	Pack of 10 x 225 ml Flask
4659	Bolton Broth Base ISO	Pack of 10 x 250 ml Flask
4675	Bolton Broth Base ISO	Pack of 10 x 90 ml Flask
5154	Buffered Peptone Water 0,1% EP/USP	Pack of 10 x 200 ml Flask
5170	Buffered Peptone Water 0,1% EP/USP	Pack of 10 x 450 ml Flask
5180	Buffered Peptone Water 0,1% EP/USP	Pack of 10 x 100 ml Flask
4638	Buffered Peptone Water 0,1% EP/USP/ISO	Pack of 10 x 90 ml Flask
4643	Buffered Peptone Water EP/USP	Pack of 10 x 90 ml Flask
4699	Buffered Peptone Water EP/USP	Pack of 10 x 225 ml Flask
5075	Buffered Peptone Water EP/USP with Neutralizings	Pack of 10 x 90 ml Flask
5020	Buffered Peptone Water ISO	Pack of 10 x 225 ml Flask
5153	Buffered Peptone Water ISO	Pack of 10 x 100 ml Flask
5171	Buffered Peptone Water ISO	Pack of 10 x 90 ml Flask
4637	Buffered Peptone water with 0,5% Polisorbate 80	Pack of 10 x 90 ml Flask
5122	Cetrimide Agar EP/USP/ISO	Pack of 10 x 100 ml Flask
5105	Chocolate Agar	Pack of 10 x 90 ml Flask
5102	CLED Agar (Cystine Lactose Electrolyte Deficient)	Pack of 10 x 100 ml Flask
5050	Columbia CNA Agar (w/o Blood)	Pack of 10 x 90 ml Flask
5077	Dey-Engley Neutralizing Broth ISO	10 x 100 ml

CAT.	PRODUCT	PACK SIZE
5044	Diluent of Beerens	Pack of 10 x 90 ml Flask
5064	Diluent of Beerens	Pack of 10 x 100 ml Flask
5065	Diluent of Beerens	Pack of 10 x 45 ml Flask
5066	Eugon LT Modified Broth	Pack of 10 x 100 ml Flask
5067	Eugon LT Modified Broth	Pack of 10 x 200 ml Flask
5068	Eugon LT Modified Broth	Pack of 10 x 90 ml Flask
5082	Fluid A EP/USP	Pack of 10 x 100 ml Flask
5083	Fluid A EP/USP	Pack of 10 X 300 ML Flask
5084	Fluid A EP/USP	Pack of 6 X 900 ML Flask
5085	Fluid D EP/USP	Pack of 10 X 100 ML Flask
5086	Fluid D EP/USP	Pack of 10 X 300 ML Flask
5087	Fluid D EP/USP	Pack of 6 X 850 ML Flask
5088	Fluid K USP	Pack of 10 X 100 ML Flask
5089	Fluid K USP	Pack of 10 X 300 ML Flask
5090	Fluid K USP	Pack of 6 X 900 ML Flask
5111	Hektoen Enteric Agar ISO	Pack of 10 x 100 ml Flask
5013	Lactose Broth	Pack of 10 x 90 ml Flask
5053	Lactose Broth	Pack of 10 x 100 ml Flask
5069	Lactose Broth	Pack of 10 x 30 ml Flask
5162	Lethen Agar Modified	Pack of 10 x 200 ml Flask
5163	Lethen Broth Modified ISO	Pack of 10 x 90 ml Flask
4697	Linden Grain Broth	Pack of 10 x 450 ml Flask
5023	Listeria Enrichment Broth Fraser ISO	Pack of 10 x 225 ml Flask
5022	Listeria Enrichment Broth Half Fraser ISO	Pack of 10 x 225 ml Flask
5003	MacConkey Agar EP/USP/ISO	Pack of 10 x 100 ml Flask
5146	Macconkey Broth EP/USP	Pack of 10 x 100 ml Flask
4692	Malt Extract Agar Acidified (pH 3,5)	Pack of 10 x 200 ml Flask
4674	Malt Extract Agar with Chloramphenicol	Pack of 10 x 200 ml Flask
5104	Mannitol Salt Agar (MSA) (Chapman Medium) EP/USP/ISO	Pack of 10 x 100 ml Flask
5093	Maximum Recovery Diluent (MRD) ISO	Pack of 10 x 250 ml Flask
4681	Milk Plate Count Agar (PCA) ISO	Pack of 10 x 200 ml Flask
5172	Milk Plate Count Agar (PCA) ISO	Pack of 10 x 100 ml Flask
5168	Mossel EE Broth EP/USP	Pack of 10 x 100 ml Flask
4691	MRS Agar	Pack of 10 x 200 ml Flask
4684	MRS Agar Low pH ISO	Pack of 10 x 200 ml Flask
5002	Mueller Hinton Agar	Pack of 10 x 100 ml Flask
4667	Muller Kauffmann Tetrathionate-Novobiocin Broth (MKTTN) ISO	Pack of 10 x 100 ml Flask
5126	Nutrient Agar ISO	Pack of 10 x 100 ml Flask
4664	OGYE Agar ISO	Pack of 10 x 100 ml Flask
4665	OGYE Agar ISO	Pack of 10 x 200 ml Flask
5113	Orange Serum Agar	Pack of 10 x 200 ml Flask
4683	Peptone Sorbitol Bile Broth (PBS) ISO	Pack of 10 x 100 ml Flask
5160	Peptone Water (Tryptone Water) ISO	Pack of 10 x 100 ml Flask
4668	Phosphate Buffered Saline (PBS)	Pack of 10 x 100 ml Flask
4676	Preston Broth without Blood	Pack of 10 x 90 ml Flask
5187	Rappaport Vassiliadis Broth EP/USP	Pack of 10 x 90 ml Flask
5107	Sabouraud Dextrose Agar EP/USP/ISO	Pack of 10 x 100 ml Flask
5143	Sabouraud Dextrose Agar EP/USP/ISO	Pack of 10 x 200 ml Flask
5108	Sabouraud Dextrose Agar with Chloramphenicol EP/USP/ISO	Pack of 10 x 100 ml Flask
4656	Sabouraud Dextrose Broth EP/USP	Pack of 10 x 100 ml Flask
5182	Saline Peptone Water ISO	Pack of 10 x 90 ml Flask

CAT.	PRODUCT	PACK SIZE
5120	Salmonella Shigella Agar (SS Agar)	Pack of 10 x 100 ml Flask
5072	Sboursaud Dextrose Agar with Chloramphenicol EP/USP/ISO	Pack of 10 x 200 ml Flask
5078	SGQ+ Broth	Pack of 10 x 450 ml Flask
5112	Standard Methods Agar (PCA) ISO/APHA	Pack of 10 x 200 ml Flask
5115	Standard Methods Agar (PCA) ISO/APHA	Pack of 10 x 100 ml Flask
4660	T.S.C. Agar Base (Tryptose Sulfite Cycloserine) ISO	Pack of 10 x 100 ml Flask
4661	T.S.C. Agar Base (Tryptose Sulfite Cycloserine) ISO	Pack of 10 x 200 ml Flask
5148	TBX Chromogenic Agar (Tryptone Bile X-Glucuronide) ISO	Pack of 10 x 100 ml Flask
5149	TBX Chromogenic Agar (Tryptone Bile X-Glucuronide) ISO	Pack of 10 x 200 ml Flask
4662	Thioglycollate Fluid Medium EP/USP/ISO	Pack of 10 x 300 ml Flask
4663	Thioglycollate Fluid Medium EP/USP/ISO	Pack of 10 x 450 ml Flask
5128	Thioglycollate Fluid Medium EP/USP/ISO	Pack of 10 x 100 ml Flask
5183	Thioglycollate Fluid Medium EP/USP/ISO	Pack of 10 x 200 ml Flask
5000	Trypticasein Soy Agar (TSA) EP/USP/ISO	Pack of 10 x 100 ml Flask
5157	Trypticasein Soy Agar (TSA) EP/USP/ISO	Pack of 10 x 200 ml Flask
4657	Trypticasein Soy Broth (TSB) EP/USP/ISO	Pack of 10 x 90 ml Flask
5070	Trypticasein Soy Broth (TSB) EP/USP/ISO	Pack of 10 x 200 ml Flask
5098	Trypticasein Soy Broth (TSB) EP/USP/ISO	Pack of 10 x 450 ml Flask
5119	Trypticasein Soy Broth (TSB) EP/USP/ISO	Pack of 10 x 100 ml Flask
4670	Violet Red Bile Agar With Glucose (VRBG) EP/USP/ISO	Pack of 10 x 200 ml Flask
5158	Violet Red Bile Agar With Glucose (VRBG) EP/USP/ISO	Pack of 10 x 100 ml Flask
4682	Violet Red Bile Lactose Agar (VRBL) ISO	Pack of 10 x 200 ml Flask
5161	Violet Red Bile Lactose Agar (VRBL) ISO	Pack of 10 x 100 ml Flask
4666	WL Nutrient Modified Agar	Pack of 10 x 100 ml Flask
4671	Yeast Extract Agar (YEA) ISO	Pack of 10 x 200 ml Flask
4672	Yeast Extract Agar (YEA) ISO	Pack of 10 x 200 ml Flask
4679	Yeast Glucose Chloramphenicol Agar (YGC) ISO	Pack of 10 x 100 ml Flask
4680	Yeast Glucose Chloramphenicol Agar (YGC) ISO	Pack of 10 x 200 ml Flask
4693	YSG Medium	Pack of 10 x 200 ml Flask

Bags

CAT.	PRODUCT	PACK SIZE
6702	Buffered Peptone Water ISO	2x5L
6705	Buffered Peptone Water ISO	3x3L
6707	Buffered Peptone Water ISO	5x2L
6703	Listeria Enrichment Broth Half Fraser ISO	2x5L
6706	Listeria Enrichment Broth Half Fraser ISO	3x3L
6708	Listeria Enrichment Broth Half Fraser ISO	5x2L
6710	Maximun Recovery Diluent (MRD)	3x3L
6711	Maximun Recovery Diluent (MRD)	2x5L
6701	Trypticasein Soy Broth (TSB) EP/ USP/ISO	2x5L
6704	Trypticasein Soy Broth (TSB) EP/USP/ISO	3x3L

Water plates 55 mm

CAT.	PRODUCT	PACK SIZE
4743	Baird Parker Agar	Pack of 30 Water plates
4701	Bile Esculin Azide Agar ISO	Pack of 30 Water plates
4739	Brettanomyces Agar	Pack of 30 Water plates
4704	Cetrimide Agar EP/USP/ISO	Pack of 30 Water plates
4702	Chapman TTC Agar (Tergitol 7 Agar)	Pack of 30 Water plates
4721	Chromogenic Coliforms Agar (CCA) ISO	Pack of 30 Water plates
4709	Clostridium Perfringens Agar (m-CP)	Pack of 30 Water plates
4729	E. coli-Coliforms Chromogenic Agar Base (BOE)	Pack of 30 Water plates
4720	Endo LES Agar	Pack of 30 Water plates
4722	Fecal Coliforms Agar (m-FC)	Pack of 30 Water plates
4550	Malt Extract Agar Acidified (pH 4,3)	Pack of 30 Water plates
4716	Mannitol Salt Agar (MSA) (Chapman Medium) EP/USP/ISO	Pack of 30 Water plates
4740	Pseudomonas Agar CN ISO	Pack of 30 Water plates
4711	R2A Agar EP/USP	Pack of 30 Water plates
4714	Sabouraud Dextrose Agar EP/USP/ISO	Pack of 30 Water plates
4713	Sabouraud Dextrose Agar with Chloramphenicol EP/USP/ISO	Pack of 30 Water plates
4710	Slanetz-Bartley Medium ISO	Pack of 30 Water plates
4712	SPS Agar (Sulfite Polymyxin Sulfadiazine)	Pack of 30 Water plates
4706	Standard Methods Agar (PCA) ISO/APHA	Pack of 30 Water plates
4728	T.S.C. Agar Base (Tryptose Sulfite Cycloserine) ISO	Pack of 30 Water plates
4703	Yeast Extract Agar (YEA) ISO	Pack of 30 Water plates

Contact plates 60 mm (Rodac)

CAT.	PRODUCT	PACK SIZE
4530	Baird-Parker Agar ISO	Pack of 30 Rodac plates
4527	Cetrimide Agar	Pack of 30 Rodac plates
4546	Dextrose Sabouraud Agar with Neutralings. Triple Wrapped. Irradiated	Pack of 24 Rodac plates
4523	MacConkey Agar EP/USP/ISO	Pack of 30 Rodac plates
4526	Mannitol Salt Agar (MSA) (Chapman Medium) EP/USP/ISO	Pack of 30 Rodac plates
4529	Rose Bengal Agar with Chloramphenicol	Pack of 30 Rodac plates
4522	Sabouraud Dextrose Agar EP/USP/ISO	Pack of 30 Rodac plates
4528	Sabouraud Dextrose Agar with Chloramphenicol	Pack of 30 Rodac plates
4544	Sabouraud Dextrose Agar with Chloramphenicol EP/USP/ISO. Trple Wrapped. Irradiated ISO	Pack of 24 Rodac plates
4545	Sabouraud Dextrose Agar with Neutralizings. Double Wrapped. Irradiated.	Pack of 30 Rodac plates
4542	Sabouraud Dextrose Agar. Triple Wrapped. Irradiated EP/USP/ISO	Pack of 24 Rodac plates
4521	Standard Methods Agar (PCA)	Pack of 30 Rodac plates
4552	Standard Methods Agar (PCA) TLHTh. Double wrap. Irradiated.	Pack of 30 Rodac plates
4553	Standard Methods Agar (PCA) TLHTh. Triple wrap. Irradiated.	Pack of 24 Rodac plates
4520	Trypticasein Soy Agar (TSA) EP/USP/ISO	Pack of 30 Rodac plates
4533	Trypticasein Soy Agar (TSA) EP/USP/ISO. Triple Wrapped. Irradiated	Pack of 24 Rodac plates
4538	Trypticasein Soy Agar (TSA) with Neutralizings. Double Wrapped. Irradiated	Pack of 30 Rodac plates
4534	Trypticasein Soy Agar (TSA) with Neutralizings. Triple Wrapped. Irradiated	Pack of 24 Rodac plates
4525	TSA Lethen Agar (Irradiated)	Pack of 30 Rodac plates
4561	TSA with Penase with Lecithin with Polisorbate	Pack of 30 Rodac plates
4524	Violet Red Bile Agar With Glucose (VRBG) EP/USP/ISO	Pack of 30 Rodac plates
4532	Violet Red Bile Lactose Agar (VRBL) ISO	Pack of 30 Rodac plates
4554	VRBG TLHTh Agar	Pack of 30 Rodac plates

Lock&Block Pharma Plates

CAT.	PRODUCT	PACK SIZE
4575	Lock&Block TSA (Trypticasein Soy Agar) TLHTh. Triple Wrapped. Irradiated EP/USP/ISO	Pack of 30 Rodac plates
4576	Lock&Block Sabouraud TLHTh Agar. Triple Wrapped. Irradiated EP/USP	Pack of 30 Rodac plates
4577	Lock&Block TSA + Penase Triple Wrapped. Irradiated	Pack of 30 Rodac plates
4578	Lock&Block Sabouraud Chloramphenicol Dextrose Agar. Triple Wrapped. Irradiated EP/USP/ISO	Pack of 30 Rodac plates
803	LC Sabouraud Dextrose Agar. Double Wrapped. Irradiated EP/USP	Pack of 20 plates
804	LC Trypticasein Soy Agar (TSA). Double wrapped. Irradiated EP/USP/ISO	Pack of 20 plates
805	LC Sabouraud TLHTh Agar. Triple Wrapped. Irradiated EP/USP	Pack of 20 plates
806	LC Sabouraud Chloramphenicol Dextrose Agar. Triple Wrapped. Irradiated EP/USP/ISO	Pack of 20 plates
807	LC Sabouraud Dextrose Agar. Triple Wrapped. Irradiated EP/USP	Pack of 20 plates
808	LC TSA + Penase Triple Wrapped. Irradiated	Pack of 20 plates
809	LC TSA (Trypticasein Soy Agar) TLHTh. Triple Wrapped. Irradiated EP/USP/ISO	Pack of 20 plates
810	LC Trypticasein Soy Agar (TSA). Triple Wrapped. Irradiated EP/USP/ISO	Pack of 20 plates

Hemocultures

CAT.	PRODUCT	PACK SIZE
3004	Brain Heart Infusion Broth with SPS, CO2 & Vacuum	Pack of 10 x 50 ml Flask
3005	Brain Heart Infusion Broth with SPS, CO2 & Vacuum	Pack 8 x 20 ml
3107	Schaedler Broth with SPS, Vacuum and CO2	Pack of 10 x 50 ml Flask
3105	Thioglycollate Broth with SPS, CO2 & Vacuum	Pack of 10 x 50 ml Flask
3001	Trypticasein Soy Broth with SPS, CO2 & Vacuum	Pack of 10 x 50 ml Flask
3008	Trypticasein Soy Broth with SPS, CO2 & Vacuum	Pack 8 x 20 ml

Mycobacteria testing

CAT.	PRODUCT	PACK SIZE
5017	B Reagent	10 Units
4100	Lowenstein Jensen Medium	Pack of 20 Tubes
4125	Lowenstein Jensen Medium (Long Tube)	Pack of 20 Tubes
4120	Lowenstein Jensen Medium with Pyruvate	Pack of 20 Tubes
4126	Lowenstein Jensen Medium with Pyruvate (Long Tube)	Pack of 20 Tubes
5015	Reagent N° 2	10 Units
6006	Reagent N°1	10 Units

CondaChrome®

CondaChrome® culture media have in their composition a colorless chromogenic substrate, which thanks to the specific enzymatic activity of each microorganism, is degraded releasing a part of it, called chromophore, giving the colony an intense and specific color that makes it possible to identify the bacteria at first sight. Discover CondaChrome®, Condalab's line of chromogenic culture media. Step into the next generation of culture media and improve the accuracy and speed of your results.

CAT.	PRODUCT	PACK SIZE
2207	Candida Chromogenic Agar	500 grams
2142	Chromogenic Agar Burkholderia Cepacia	500 grams
2080	Chromogenic Coliforms Agar (CCA) ISO	500 grams
1446	Chromogenic Cronobacter Isolation Agar (CCI) ISO	500 grams
1588	E. coli O157:H7 Chromogenic Agar Base	500 grams
1491	E. coli-Coliforms Chromogenic Agar Base (BOE)	500 grams
1340	E. coli-Coliforms Chromogenic Medium	500 grams
2018	E. coli-Enterobacteria Chromogenic Agar	500 grams
1285	EC with MUG Fluorogenic Broth	500 grams
2062	ESBL Chromogenic Agar	500 grams
2119	Klebsiella Chromogenic Agar Base	500 grams
2063	KPC Chromogenic Medium	500 grams
1465	Lauryl Sulfate Chromogenic Broth	500 grams
1345	Listeria Chromogenic Agar Base according to Ottaviani and Agosti (ALOA) ISO	500 grams
1412	M-EI Chromogenic Agar Base	500 grams
2050	m-EI Chromogenic Agar Base, Modified	500 grams
1423	MRSA Chromogenic Agar Base	500 grams
1498	MRSA Chromogenic Modified Agar Base	500 grams
2144	PEC Chromogenic Agar	500 grams
1493	Pseudomonas Chromogenic Agar	500 grams
1122	Salmonella Chromogenic Agar	500 grams
1585	Standard Method Chromogenic Agar (PCA)	500 grams
2076	Staphylococcus Chromogenic Agar	500 grams
1151	TBX Chromogenic Agar (Tryptone Bile X-Glucuronide) ISO	500 grams
1424	Urinary Tract Infections Chromogenic Agar (UTIC)	500 grams
2077	Vancomycin-Resistant Enterococcus (VRE) Chromogenic Agar	500 grams
2054	Vibrio Chromogenic Agar	500 grams

Supplements and additives.

In order to increase a given culture medium selectivity, a supplement or additive not included in the medium (due to its sensitivity to sterilization temperature) must be added on certain occasions, to ensure the isolation and identification of specific microorganisms. Condalab supplements are now available in freeze-dried vials and in vials with sterile solutions, without the need for previous dissolution.

CAT.	PRODUCT	PACK SIZE
6052	Ampicillin Supplement (Aeromonas)	Pack of 10 vials for 500 ml/each
6021	Bacillus Cereus Supplement ISO	Pack of 10 vials for 500 ml/each
6070	Bolton Selective Supplement	Pack of 10 vials for 500 ml/each
6015	Bordetella Supplement	Pack of 10 vials for 500 ml/each
6017	Brucella Supplement	Pack of 10 vials for 500 ml/each
6060	Brucella Supplement	Pack of 10 vials for 500 ml/each
6032	Burkholderia cepacia Supplement	Pack of 10 vials for 500 ml/each
6053	CCDA Supplement (Campylobacter Blood free)	Pack of 10 vials for 500 ml/each
6064	Cefixime Tellurite (CT) Supplement	Pack of 10 vials for 500 ml/each
6069	Cefoxitin Supplement	Pack of 10 vials for 500 ml/each
6036	CFC Supplement	Pack of 10 vials for 500 ml/each
6049	Chloramphenicol Selec. Suppl. (25 Mg)	Pack of 10 vials for 500 ml/each
6027	Chloramphenicol Supplement	Pack of 10 vials for 500 ml/each
6061	Clostridium Difficile Supplement	Pack of 10 vials for 500 ml/each
6020	Clostridium Perfringens Supplement (TSC)	Pack of 10 vials for 500 ml/each
6016	CNA Supplement (Staph/Sterp)	Pack of 10 vials for 500 ml/each
6041	E. coli-Coliforms Supplement (BOE)	Pack of 10 vials for 500 ml/each
5152	Egg Yolk Emulsion	100 ml Flask
5156	Emulsión Yema de Huevo 50%	500 ml
6042	ESBL Supplement	Pack of 10 vials for 500 ml/each
6023	Fecal Coliforms Supplement	Pack of 10 vials for 500 ml/each
6050	Ferric Ammonium Citrate Supplement	Pack of 10 vials for 500 ml/each
6001	Fraser Listeria Selective Supplement	Pack of 10 vials for 500 ml/each
6002	Half Fraser Listeria Selective Supplement	Pack of 10 vials for 500 ml/each
6051	ITC Supplement	Pack of 10 vials for 480 ml/each
6045	Klebsiella Selective Supplement	Pack of 10 vials for 500 ml/each
5205	Kovac's Reagent	100 ml Flask
6012	LCAT Supplement	Pack of 10 vials for 500 ml/each
6022	Legionella BCYE Supplement ISO	Pack of 10 vials for 500 ml/each
6092	Legionella BCYE-cys supplement ISO	2x5 vials
6025	Legionella GVPC Supplement ISO	Pack of 10 vials for 500 ml/each
6067	Legionella MWY	Pack of 10 vials for 500 ml/each
6040	Listeria Chromogenic Selective Supplement	Pack of 10 vials for 500 ml/each
6031	Listeria Lipase C Supplement	Pack of 10 vials for 500 ml/each
6073	m-CP Supplement	Pack of 10 vials for 500 ml/each
6046	Nalidixic Selective Supplement	Pack of 10 vials for 500 ml/each
6048	Novobiocin Select. Suppl. (10 Mg)	Pack of 10 vials for 500 ml/each
6047	Novobiocin Select. Suppl. (20 Mg)	Pack of 10 vials for 1000 ml/each
6018	OGA Supplement	Pack of 10 vials for 500 ml/each
6003	Oxford Listeria Selective Supplement	Pack of 10 vials
6007	Oxidase Reagent	2 x 25 swabs
6004	Palcam Listeria Selective Supplement	Pack of 10 vials
6011	Polyenrichment Supplement	Pack of 2x5 vials for 500 ml/pair vials
6090	Polyenrichment Supplement (Each pair of vial for 5L of medium)	2x4 vials

CAT.	PRODUCT	PACK SIZE
5209	Potassium Tellurite 1%	100 ml Flask
5208	Potassium Tellurite 3,5%	100 ml Flask
6019	Preston Campylobacter Supplement	Pack of 10 vials for 500 ml/each
6081	Preston Campylobacter Supplement ISO	Pack of 10 vials for 500 ml/each
6024	RPF Supplement ISO	Pack of 10 vials for 100 ml/each
6043	Salmonella Chromogenic Supplement	Pack of 10 vials for 500 ml/each
6089	Selective Supplement for Burkholderia cepacia (BCSA)	Pack of 10 vials for 500 ml/each
6074	Selective Supplement MUP	Pack of 10 vials for 500 ml/each
5129	Tellurite Egg Yolk Emulsion	100 ml Flask
5165	Tellurite Egg Yolk Emulsion 250 ml	250 ml Flask
5164	Tellurite Egg Yolk Emulsion 500 ml	500 ml Flask
6034	Tinsdale Enrichment Supplement	Pack of 20 tubes for 15 ml/each
6030	TTC 1% Supplement	Pack 10 vials (1 vial / 500 ml of medium)
5100	Urea 40% Solution	100 ml Flask
6072	Vancomycin Supplement	Pack of 10 vials for 10 ml/each
6014	VCAT Supplement	Pack of 10 vials for 500 ml/each
6013	VCN Supplement	Pack of 10 vials for 500 ml/each
6026	VCNT Supplement	Pack of 10 vials for 500 ml/each
6062	XLT4 Supplement	100 ml Flask
6033	Yersinia Selective Supplement (CIN)	Pack of 10 vials for 500 ml/each

Bacteriological stains.

Stains and dyes are mainly used in procedures for microorganism differentiation and identification, based on specific characteristics of bacteria. This step is of vital importance when identifying a microorganism and establishing a reliable clinical diagnosis in microbiology and pathological anatomy.

CAT.	PRODUCT	PACK SIZE	CAT.	PRODUCT	PACK SIZE
4614	Acetone Ethanol Decolorant	1 x 1L	4603	Methylene Blue	1 x 1L
5059	Acetone Ethanol Decolorant	1 x 250 ml	5058	Methylene Blue	1 x 250 ml
4647	Auramide	1 x 1L	4610	Oxalate Crystal Violet	1 x 1L
5005	Auramide	1 x 250 ml	5062	Oxalate Crystal Violet	1 x 250 ml
4611	Chlorhidric Alcohol	1 x 1L	4655	Permanganate	1 x 1L
5024	Chlorhidric Alcohol	1 x 250 ml	5014	Permanganate	1 x 250 ml
4618	Fenicade Fuchsine	1 x 1L	4616	Safranin	1 x 1L
5060	Fenicade Fuchsine	1 x 250 ml	5063	Safranin	1 x 250 ml
4600	Gram Stain KIT	Pack of 4 x 250 ml	4632	Thiazine Red	1 x 250 ml
4612	Lugol	1 x 1L	4651	Thiazine Red	1 x 1L
5061	Lugol	1 x 250 ml			

Antimicrobial susceptibility testing.

The analysis of microorganism sensitivity to different antimicrobial agents is a key point in any clinical microbiology laboratory. This analysis is carried out by means of an antibiogram; the objective of which is to evaluate in a laboratory the response of a microorganism to battery of antibiotics or antifungals.

Antibiotic discs.

These discs (6 mm) are prepared by impregnating high quality absorbent paper with accurately determined amounts of antibiotic or other chemotherapeutic agents. It is used in antibiogram techniques which allows to in vitro subject a standard bacterial inoculum to one or different antibiotic concentration to predict in vivo effectiveness.

CAT. N°	PRODUCT	PACK SIZE
7036	Amikacin 10 µg	5x50 Discs
7000	Amikacin 30 µg	5x50 Discs
7086	Amoxicillin / Clavulanic acid 3 µg	5x50 Discs
7002	Amoxicillin / Clavulanic acid 30 µg	5x50 Discs
7038	Amoxicillin 10 µg	5x50 Discs
7001	Amoxicillin 25 µg	5x50 Discs
7085	Amoxicillin 30 µg	5x50 Discs
7087	Amphotericin B 20 U	5x50 Discs
7005	Ampicillin / Sulbactam 20 µg	5x50 Discs
7152	Ampicillin / Sulbactam 30 µg	5x50 Discs
7004	Ampicillin 10 µg	5 x 50 disc
7137	Ampicillin 2 µg	5x50 Discs
7142	Ampicillin 25 µg	5x50 Discs
7095	Ampicillin 30 µg	5x50 Discs
7006	Azithromycin 15 µg	5x50 Discs
7007	Aztreonam 30 µg	5x50 Discs
7008	Bacitracin 0.04 µg	5x50 Discs
7157	Bacitracin 10 U	5x50 Discs
7387	Blank Disc	5x50 Discs
7097	Carbenicillin 100 µg	5x50 Discs
7098	Cefadroxil 30 µg	5x50 Discs
7010	Cefamandole 30 µg	5x50 Discs
7011	Cefazolin 30 µg	5x50 Discs
7013	Cefepime / Clavulanic acid 40 µg	5x50 Discs
7012	Cefepime 30 µg	5x50 Discs
7168	Cefixime 10 µg	5x50 Discs
7015	Cefixime 5 µg	5x50 Discs
7171	Cefoperazone / Sulbactam 105 µg	5x50 Discs
7170	Cefoperazone 30 µg	5x50 Discs
7016	Cefoperazone 75 µg	5x50 Discs
7020	Cefotaxime / Clavulanic acid 40 µg	5x50 Discs
7101	Cefotaxime 10 µg	5x50 Discs
7017	Cefotaxime 30 µg	5x50 Discs
7018	Cefotaxime 5 µg	5x50 Discs
7022	Cefoxitin 30 µg	5x50 Discs
7436	Ceftazidime / Avibactam 14 µg	5x50 Discs
7027	Ceftazidime / Clavulanic acid 40 µg	5x50 Discs
7103	Ceftazidime 10 µg	5x50 Discs
7026	Ceftazidime 30 µg	5x50 Discs

CAT. N°	PRODUCT	PACK SIZE
7028	Ceftizoxime 30 µg	5x50 Discs
7029	Ceftriaxone 30 µg	5x50 Discs
7186	Ceftriaxone 5 µg	5x50 Discs
7030	Cefuroxime 30 µg	5x50 Discs
7031	Cephalexin 30 µg	5x50 Discs
7032	Cephalothin 30 µg	5x50 Discs
7033	Chloramphenicol 30 µg	5x50 Discs
7199	Ciprofloxacin 10 µg	5x50 Discs
7200	Ciprofloxacin 30 µg	5x50 Discs
7034	Ciprofloxacin 5 µg	5x50 Discs
7107	Clindamycin 10 µg	5x50 Discs
7040	Clindamycin 2 µg	5x50 Discs
7041	Clotrimazole 10 µg	5x50 Discs
7210	Cloxacillin 5 µg	5x50 Discs
7043	Colistin 10 µg	5x50 Discs
7213	Colistin 25 µg	5x50 Discs
7214	Colistin 50 µg	5x50 Discs
7220	Doxycycline 10 µg	5x50 Discs
7044	Doxyxyciclina 30 µg	5x50 Discs
7221	Econazole 10 µg	5x50 Discs
7233	Enrofloxacin 5 µg	5x50 Discs
7047	Ertapenem 10 µg	5x50 Discs
7048	Erythromycin 15 µg	5x50 Discs
7235	Erythromycin 30 µg	5x50 Discs
7238	Florfenicol 30 µg	5x50 Discs
7241	Fluconazole 25 µg	5x50 Discs
7049	Fosfomicin 200 µg	5x50 Discs
7114	Fosfomicin 50 µg	5x50 Discs
7247	Furazolidone 100 µg	5x50 Discs
7249	Furazolidone 50 µg	5x50 Discs
7050	Fusidic acid 10 µg	5x50 Discs
7051	Gentamicin 10 µg	5x50 Discs
7117	Gentamicin 120 µg	5x50 Discs
7118	Gentamicin 30 µg	5x50 Discs
7052	Imipenem 10 µg	5x50 Discs
7120	Itraconazole 10 µg	5x50 Discs
7053	Kanamycin 30 µg	5x50 Discs
7269	Ketoconazole 10 µg	5x50 Discs
7271	Levofloxacin 15 µg	5x50 Discs

CAT. N°	PRODUCT	PACK SIZE
7056	Levofloxacin 5 µg	5x50 Discs
7057	Lincomycin 15 µg	5x50 Discs
7273	Lincomycin 2 µg	5x50 Discs
7122	Linezolid 10 µg	5x50 Discs
7121	Linezolid 30 µg	5x50 Discs
7279	Marbofloxacin 5 µg	5x50 Discs
7059	Meropenem 10 µg	5x50 Discs
7124	Metronidazole 5 µg	5x50 Discs
7060	Miconazole 10 µg	5x50 Discs
7062	Moxifloxacin 5 µg	5x50 Discs
7128	Mupirocin 5 µg	5x50 Discs
7063	Nalidixic acid 30 µg	5x50 Discs
7383	Neomycin 10 µg	5x50 Discs
7064	Neomycin 30 µg	5x50 Discs
7296	Netilmicin 10 µg	5x50 Discs
7065	Netilmicin 30 µg	5x50 Discs
7131	Nitrofurantoin 100 µg	5x50 Discs
7066	Nitrofurantoin 300 µg	5x50 Discs
7301	Nitroxoline 30 µg	5x50 Discs
7067	Norfloxacin 10 µg	5x50 Discs
7305	Novobiocin 30 µg	5x50 Discs
7068	Novobiocin 5 µg	5x50 Discs
7069	Nystatin 100 U µg	5x50 Discs
7070	Ofloxacin 10 µg	5x50 Discs
7071	Ofloxacin 5 µg	5x50 Discs
7073	Optochin 5 µg	5x50 Discs
7074	Oxacillin 1 µg	5x50 Discs
7313	Oxacillin 5 µg	5x50 Discs
7315	Oxytetracycline 30 µg	5x50 Discs
7076	Pefloxacin 5 µg	5x50 Discs
7134	Penicillin G 1 U	5x50 Discs
7077	Penicillin G 10 U	5x50 Discs
7320	Phenoxyethylpenicillin 10 µg	5x50 Discs
7079	Piperacillin / Tazobactam 110 µg	5x50 Discs
7136	Piperacillin / Tazobactam 36 µg	5x50 Discs
7326	Piperacillin / Tazobactam 85 µg	5x50 Discs
7078	Piperacillin 100 µg	5x50 Discs
7135	Piperacillin 30 µg	5x50 Discs
7080	Polymixin B 300 U µg	5x50 Discs
7081	Rifampin (Rifampicin) 5 µg	5x50 Discs
7082	Spectinomycin 100 µg	5x50 Discs
7347	Spiramycin 100 µg	5x50 Discs
7083	Streptomycin 10 µg	5x50 Discs
7349	Streptomycin 25 µg	5x50 Discs
7084	Streptomycin 300 µg	5x50 Discs
7360	Sulphamethoxazole 50 µg	5x50 Discs
7088	Teicoplanin 30 µg	5x50 Discs
7089	Telithromycin 15 µg	5x50 Discs
7362	Tetracycline 10 µg	5x50 Discs
7091	Tetracycline 30 µg	5x50 Discs

CAT. N°	PRODUCT	PACK SIZE
7144	Ticarcillin / Clavulanic acid 85 µg	5x50 Discs
7092	Ticarcillin 75 µg	5x50 Discs
7093	Tigecycline 15 µg	5x50 Discs
7094	Tobramycin 10 µg	5x50 Discs
7365	Tobramycin 30 µg	5x50 Discs
7352	Trimethoprim / Sulfadiazine 30 µg	5x50 Discs
7042	Trimethoprim / Sulphamethoxazol 25 µg	5x50 Discs
7366	Trimethoprim 1,25 µg	5x50 Discs
7370	Trimethoprim 5 µg	5x50 Discs
7376	Tylosin 30 µg	5x50 Discs
7378	V Factor Disc	5x50 Discs
7145	Vancomycin 30 µg	5x50 Discs
7380	Vancomycin 5 µg	5x50 Discs
7159	X + V factor disc	5x50 Discs
7160	X Factor Disc	5x50 Discs

Ingredients.

Over 60 years of experience makes Condalab a true expert in the whole development process of culture media which is why we also offer all the necessary raw materials for the manufacturing of media.

Agars.

Agar is a polysaccharide that can be found in some marine algae cell walls. Condalab's agar is obtained from those species that provide a higher quality performance such as Gelidium algae. Agar is the gelling agent in solid culture media. The high-quality agar offered by Condalab has no impurities and is microbiologically neutral, to obtain the best results in culture media preparation and other applications.

CAT.	PRODUCT	PACK SIZE	CAT.	PRODUCT	PACK SIZE
1802	American Bacteriological Agar	500 grams	1816	Pharmaceutical Agar	500 grams
1800	European Bacteriological Agar	500 grams	1812	Plant Propagation Agar	500 grams
1804	Industrial Agar	500 grams	1806	Purified Agar	500 grams

Peptones.


Peptones, protein of vegetable or animal origin digests, constitute a mixture of amino acids and polypeptides easily assimilated by microorganisms. In its catalog, Condalab offers a wide range of peptones and extracts to be used in the manufacturing of microbiological culture media and bacterial fermentation's.

CAT.	PRODUCT	PACK SIZE	CAT.	PRODUCT	PACK SIZE
1604	Acid Casein Peptone	500 grams	1600	Meat Peptone	500 grams
1704	Bacteriological Gelatin	500 grams	1628	Peptonized Milk	500 grams
1710	Bacteriological Ox Bile	500 grams	1610	Polypeptone	500 grams
1616	Bacteriological Peptone	500 grams	1712	Pork Brain Heart Infusion	500 grams
1700	Beef Extract	500 grams	1716	Pork Heart Infusion	500 grams
1706	Bile Salts N°3	500 grams	1624	Pork Meat Peptone	500 grams
1714	Bovine Heart Infusion	500 grams	1609	Proteose peptone	500 grams
1603	Casein CC Peptone	500 grams	1607	Proteose Peptone N°3	500 grams
1602	Casein Peptone	500 grams	1615	Soy Peptone GMO-Free and Animal-Free	500 grams
1606	Gelatin Peptone	500 grams	1612	Tryptone	500 grams
1935	Hemoglobin	500 grams	1614	Tryptose	500 grams
1626	Lactoalbumin Hydrolysate	500 grams	1736	Yeast Extract	500 grams
1703	Liver Peptone	500 grams	1702	Yeast extract	500 grams
1708	Malt Extract	500 grams	1721	Yeast Extract Bioprocess Grade	500 grams

Carbohydrates.

Carbohydrates are the main source of carbon in culture media. They are essential to provide a high nutritional value and even act as selective agents depending on the capacity of microorganisms to ferment.

CAT.	PRODUCT	PACK SIZE
1900	Dextrose	500 grams
1902	Lactose	500 grams
1904	Maltose Certified	500 grams
1906	Sucrose	500 grams



**MOLECULAR
BIOLOGY**

Agaroses.

Agarose is a natural polysaccharide that forms an inert matrix used for biomolecule purification and protein separation, DNA and other macromolecules by electrophoresis.

Standard agaroses.

CAT.	PRODUCT	PACK SIZE
8024	Agarose D1 High EEO	50, 100, 250 and 500 g
8010	Agarose D1 Low EEO	50, 100, 250 and 500 g
8017	Agarose D1 Low EEO GQT	50, 100, 250 and 500 g
8019	Agarose D1 Medium EEO	50, 100, 250 and 500 g
8032	Agarose D2	50, 100, 250 and 500 g
8100	Agarose E	50, 100, 250 and 500 g

Agarose for DNA recovery.

CAT.	PRODUCT	PACK SIZE
8050	Agarose LM	50, 100, 250 and 500 g
8091	Agarose LM GQT	50, 100, 250 and 500 g
8092	Agarose LM Sieve	50, 100, 250 and 500 g
8093	Agarose NovaGel GQT	50, 100, 250 and 500 g

Agarose for high DNA resolution.

CAT.	PRODUCT	PACK SIZE
8067	Agarose MS - 12	50, 100, 250 and 500 g
8075	Agarose MS - 4	50, 100, 250 and 500 g
8001	Agarose MS - 6 Metagel	50, 100, 250 and 500 g
8065	Agarose MS - 8	50, 100, 250 and 500 g

Special agarose.

CAT.	PRODUCT	PACK SIZE
8090	Agarose FP DNA	50, 100, 250 and 500 g
8045	Agarose D-5	50, 100, 250 and 500 g

Agarose reagents.

CAT.	PRODUCT	PACK SIZE
4687	CondaSafe Stain	1 ml
CK130	PronaSafe Stain	1 ml

Low melting agaroses.

CAT.	PRODUCT	PACK SIZE
8050	Agarose LM	50, 100, 250 and 500 g
8092	Agarose LM Sieve	50, 100, 250 and 500 g
8093	Agarose NovaGel GQT	50, 100, 250 and 500 g

Molecular screening.

CAT.	PRODUCT	PACK SIZE
8091	Agarose LM GQT	50, 100, 250 and 500 g
8067	Agarose MS - 12	50, 100, 250 and 500 g
8075	Agarose MS - 4	50, 100, 250 and 500 g
8001	Agarose MS - 6 Metagel	50, 100, 250 and 500 g
8065	Agarose MS - 8	50, 100, 250 and 500 g

Selection antibiotics.

Selection antibiotics for the addition to specific microbiological culture media of molecular techniques such as cloning or bacterial transformation.

CAT.	PRODUCT	PACK SIZE
6800	Ampicillin (Sodium Salt)	5 grams
6801	Ampicillin (Sodium Salt)	25 grams
6802	Ampicillin (Sodium Salt)	50 grams
6803	Carbenicillin Disodium	5 grams
6804	Chloramphenicol	25 grams
6805	Gentamicin	5 grams
6806	Kanamycin (Monosulphate)	5 grams
6807	Kanamycin (Monosulphate)	25 grams

Microbiological culture media.

Classic concepts such as microbiology are still essential for the daily development of a molecular biology laboratory in a cutting-edge field of research.

CAT.	PRODUCT	PACK SIZE
1513	YP Agar Base Medium	500 g
1540	Induction Base Medium	500 g
1543	Grace's Medium Supplement	100 g
1560	Aminoacid Mix	100 g
2090	LB Agar with Kanamycin 50µg/ml (Lennox)	500 g
2091	Luria Agar with Kanamycin 50µg/ml (Miller's LB agar)	500 g
2093	2YT Autoinducible Growth Medium w/o Trace Elements	500 g
2112	Super Growth Autoinducible with Trace Elements	500 g
2113	2YT Autoinducible with Trace Elements	500 g
1167	2xYT Agar	500 g
1507	2xYT Medium	500 g
2095	AIM -Super Growth w/o Trace Element	500 g
1083	LB Agar (Lennox)	500 g
2084	LB Agar (Lennox) with Ampicilin 100 µg/ml	500 g
1231	LB Broth (Lennox)	500 g
2094	LB Broth Autoinducible w/o Trace Element	500 g
2114	LB Broth Autoinducible with Trace Element	500 g
1552	Luria Agar (Miller's LB Agar)	500 g
1308	Luria Agar (Miller's Modification)	500 g
2085	Luria Agar with Ampicillin 100 µg/ml (Miller's LB Agar)	500 g
2086	Luria Agar with Ampicillin 50 µg/ml (Miller's LB Agar)	500 g
2092	Luria Agar with Chloramphenicol 34 µg/ml (Miller's LB Agar)	500 g
1551	Luria Broth (Miller's LB Broth)	500 g
1266	Luria Broth (Miller's Modification)	500 g
1548	Modified Terrific Broth	500 g
1549	NZCYM Broth	500 g
1542	RM Agar Base	500 g
1538	RM Base Medium	500 g
1432	sLB Agar	500 g
1163	sLB Broth	500 g
1199	sLB Broth (Buffered)	500 g
1541	SOB Medium	500 g
2019	SOC Medium	500 g
1246	Terrific Broth	500 g
1499	Terrific Broth (Non animal origin)	500 g
2072	Terrific Broth Base Autoinducible w/o Metals	500 g
2111	Terrific Broth Base Autoinducible with Trace Element	500 g
2039	Vegetable LB Agar (Lennox)	500 g
2038	Vegetable LB Broth (Lennox)	500 g
1553	Yeast Nitrogen Base w/o Added Aminoacids and w/o Ammonium Sulfate	500 g
1545	Yeast Nitrogen Base w/o Amino Acids	500 g
1511	YP Base Medium	500 g
1546	YPD Agar	500 g
1547	YPD Broth	500 g



PCR - REAL TIME DETECTION KITS



The detection of pathogens by PCR-real time methods is a trend that has taken off rapidly in recent years. These techniques allow the detection of microorganisms through their genetic material and have an unbeatable sensitivity and specificity. Moreover, thanks to these methods we can detect pathogens in record time at low concentrations. Condagene® range has everything you need to perform your qPCR assays.

Take the step into the future of pathogen detection with Condagene®

Extraction kits.

CAT.	PRODUCT	PACK SIZE
6500	Extraction Kit Condagene® Complex	100 rnx
6501	Extraction Kit Condagene® Quick	100 rnx

Pathogen detection kits.

CAT.	PRODUCT	PACK SIZE
6517	Condagene® Legionella spp	100 rnx
6516	Condagene® Listeria monocytogenes	100 rnx
6515	Condagene® Salmonella spp	100 rnx



CELL CULTURE

Condalow®

Metabolites, biomolecules, and biomass manufacturing by cell culture is an extremely delicate process. Only the best ingredients and supplements can ensure an efficient and safe production. Condalow® range give to you the highest quality peptones with low endotoxin levels and higher nutritional value because, unlike other production processes, they do not undergo ultrafiltration.

Optimize the growth and production of your cell cultures thanks to Condalow®

Condalow®

CAT.	PRODUCT	PACK SIZE
1636	Condalow® Casein Peptone	25 Kg
1637	Condalow® Meat Peptone	25 Kg
1618	Condalow® Soy Peptone	25 Kg
1638	Condalow® Wheat Peptone	25 kg

Additional information.

General storage conditions ^{1,2,3}

2 °C to 25 °C

Agaroses
Agars
Bacteriological stains
Carbohydrates
Chromogenic media
Dehydrated culture media
Hemocultures
Microbiological culture media for molecular biology
Mycobacteria testing tubes
Peptones
Selection antibiotics

8 °C to 25 °C

Bags
Contact plates 60 mm (Rodac)
Flasks
Tubes
Water plates 55 mm

2 °C to 14 °C

Plates 90 mm (Petri dish)
Supplements and additives

-20 °C to 8 °C

Antibiotic discs

General expiry dates / retest data ^{1,2,3}

2.5 to 3.5 months

Plates 90 mm (Petri dish)

6 to 7 months

Antibiotic discs
Water plates 55 mm

9 to 12 months

Tubes (including Mycobacteria testing tubes)

12 to 16 months

Bags
Flasks

9 to 18 months

Bacteriological stains
Hemocultures

12 to 24 months

Chromogenic media
Selection antibiotics

24 to 48 months

Agaroses
Agars
Carbohydrates
Dehydrated culture media
Microbiological culture media for Molecular biology
Peptones
Supplements and additives

Inspired by knowledge

1. There may be exceptions in these categories

2. Please check the shelf life and storage conditions of each reference with your Area Manager/Sales Representative


3. Original sealed packaging indicates storage and shelf life



Condalab

Inspired by knowledge

export@condalab.com | www.condalab.com



www.condalab.com