



FineTest[®] ELISA Kit

your good research partner

Wuhan Fine Biotech Co., Ltd

+86-027-87384275/87800889 +86-027-87800889

fine@fn-test.com www.fn-test.com

B9 Bld, High-Tech Medical Devices Park, No.818 Gaoxin Ave.
East Lake High-Tech Development Zone, Wuhan, China(430206)



ELISA Kits Feature and Advantage



FineTest® is dedicated to promoting global life science research by providing a diverse selection with high quality and affordable pricing also for quite extensive species diversity;

Our kits can serve a broad range of research area, including Cancer, Neuroscience, Signal Transduction, Cardiovascular, Developmental biology, Cell Division and Proliferation, Immunology and Metabolism. In addition, our kits are applicable to detection of a wide variety of substances, including proteins, peptides, and small molecules.

ELISA Kits Feature and Advantage

Advantages

- Complete and ready- to- use ELISA, incorporated all reagents and components for experiment.
- High specificity: Validated by various cross-reactivity tests, avoid interference of homology.
- Precision: recovery and linearity validated to ensure the precision and validity of test results.
- Reproducibility: strict QC control for each procedure to keep CV deviation under a very low level and ensure good reproducibility.

ELISA Kit Components

- ☑ Pre-coated ELISA Microplate
- ☑ Biotin-Labeled Antibody
- ☑ HRP-Streptavidin Conjugate
- ☑ Standard
- ☑ Optimized Sample Dilution Buffer
- ☑ Wash Buffer
- ☑ TMB Substrate
- ☑ Stop Solution
- ☑ Plate Sealer
- ☑ Product Description

EH0201 Human IL-6 ELISA Kit Data

Precision

- Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision.
- Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision.

Sample ID	Intra assay precision			Intra assay precision		
	1#	2#	3#	1#	2#	3#
n	20	20	20	20	20	20
Mean(pg/ml)	210	113	61	203	105	57
Standard deviation	3.2	1.4	0.6	6.2	3.0	1.3
CV(%)	2.8	3.5	4.1	3.0	4.1	5.6

Sensitivity

Human IL-6 ELISA kit the minimum detection range can be at 0.12 pg/ml.

ELISA Kits Feature and Advantage

Nature Proteins Testing Data

Human blood mononuclear cells (1.5×10^6 cells/mL) were cultured in RPMI supplemented with 10% fetal bovine serum, 2mM/L glutamine, 100U/mL penicillin, and 100µg/mL streptomycin sulfate and stimulated for 1,3,5 days with 15µg/mL PHA-L, assayed for levels of human IL-6 protein.

	Day 1	Day 3	Day 5
Unstimulated	375pg/ml	411pg/ml	407pg/ml
Stimulated	8322pg/ml	11,632pg/ml	13,370pg/ml

Linearity

Spiked human negative serum and cell culture medium(negative) samples to recombinant human IL-6 were serially diluted and assayed for recovery, values between 92-110% show good assay linearity.

Serum samples	Dilution	Average of Expected	Range
N=10	1:2	95-110%	99%
	1:4	91-107%	95%
	1:8	88-115%	92%
Cell culture samples	Dilution	Average of Expected	Range
N=10	1:2	97-120%	100%
	1:4	92-117%	97%
	1:8	90-119%	93%

Specificity

Test the cross relation of GM-CSF, IL-1α, IL-1β, IL-2, IL-3, IL-4, IL-8(Standard protein concentration 20ng/ml) between human and mouse. There is not any interference or positive signal.

Stability

The stability of ELISA Kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% within the expiration date under appropriate storage condition.

Standard(n=5)	37°C for 1 month	4°C for 6 months
Average (%)	80	95-100

Hot Selling Kits

Cat. No.	Product Name	Sensitivity	Detection Range
EH0024	Human VEGFR1(Vascular Endothelial Growth Factor Receptor 1) ELISA Kit	75pg/ml	125-8000pg/ml
EH0161	Human sICAM-1/CD54(Soluble Intercellular Adhesion Molecule 1) ELISA Kit	0.469ng/ml	0.781-50ng/ml
EH0209	Human SERPINA4(Kallistatin) ELISA Kit	9.375ng/ml	15.625-1000ng/ml
EH0124	Human SELE(E-Selectin) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EH0075	Human sCD163(Soluble Cluster of Differentiation 163) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EH0062	Human PRB4(Basic Salivary Proline Rich Protein 4) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EH0055	Human PRB3(Basic Salivary Proline Rich Protein 3) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0022	Human PLGF(Placental Growth Factor) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0206	Human PIM-1(Serine/Threonine-protein Kinase Pim-1) ELISA Kit	0.180ng/ml	0.313-20ng/ml
EH0021	Human PECAM-1(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit	0.469ng/ml	0.781-50ng/ml
EH0242	Human NGF/NGFβ(Beta-nerve growth factor) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
EH0018	Human MPO(Myeloperoxidase) ELISA Kit	0.938ng/ml	1.56-100ng/ml
EH0238	Human MMP-9(Matrix Metalloproteinase 9) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH0237	Human MMP-8(Matrix Metalloproteinase 8) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0017	Human MMP-2(Matrix Metalloproteinase 2) ELISA Kit	0.191ng/ml	0.5-32ng/ml
EH0234	Human MMP13(Matrix Metalloproteinase 13) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0232	Human MMP-1(Matrix Metalloproteinase 1) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EH0070	Human MME(Nepirlysin) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0002	Human MIP-3β(Macrophage Inflammatory Protein 3β) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0015	Human MICA(MHC class I polypeptide-related sequence A) ELISA Kit	46.9pg/ml	78-5000pg/ml
EH0223	Human MDC(Macrophage Derived Chemokine) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EH0014	Human M-CSF(Macrophage Colony-Stimulating Factor) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EH0222	Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0214	Human LN(Laminin) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0012	Human Lipocalin-2(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH0216	Human LEP(Leptin) ELISA Kit	18.75pg/ml	31.2-2000pg/ml
EH0208	Human KLK6(Kallikrein-6) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH0210	Human KIM-1(Kidney Injury Molecule 1) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0205	Human IL-8(Interleukin 8) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EH0201	Human IL-6(Interleukin 6) ELISA Kit	2.813pg/ml	4.69-300pg/ml
EH0200	Human IL-5(Interleukin 5) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0199	Human IL-4(Interleukin 4) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EH0197	Human IL-31(Interleukin-31) ELISA Kit	4.688pg/ml	7.813-500pg/ml
EH0196	Human IL-3(Interleukin 3) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0193	Human IL-27(Interleukin-27 subunit alpha) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0191	Human IL-22(Interleukin-22) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0189	Human IL-2(Interleukin 2) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
EH0188	Human IL1RL1(Interleukin-1 receptor-like 1) ELISA Kit	46.9pg/ml	78-5000pg/ml
EH0187	Human IL1R2(Interleukin 1 Receptor Type II) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH0185	Human IL-1beta(Interleukin 1 Beta) ELISA Kit	2.344pg/ml	3.9-250pg/ml
EH0011	Human IL-18(Interleukin 18) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0180	Human IL-17E/IL-25(Interleukin 25) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
EH0177	Human IL-15(Interleukin 15) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0176	Human IL-12(Interleukin 12) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0174	Human IL-11(Interleukin 11) ELISA Kit	18.75pg/ml	31.25-2000pg/ml

Hot Selling Kits

Cat. No.	Product Name	Sensitivity	Detection Range
EH0173	Human IL-10(Interleukin 10) ELISA Kit	4.688pg/ml	7.813-500pg/ml
EH0171	Human IGFBP-7(Insulin-like growth factor-binding protein 7) ELISA Kit	0.563ng/ml	0.938-60ng/ml
EH0170	Human IGFBP-6(Insulin-like growth factor-binding protein 6) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0169	Human IGFBP-3(Insulin-like growth factor-binding protein 3) ELISA Kit	0.234ng/ml	0.391-25ng/ml
EH0166	Human IGF-2(Insulin-like growth factor II) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EH0165	Human IGF-1(Insulin-Like Growth Factor 1) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EH0164	Human IFN γ (Interferon gamma) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0084	Human IFN κ (Interferon Kappa) ELISA Kit	4.69pg/ml	7.8-500pg/ml
EH0159	Human HGF(Hepatocyte Growth Factor) ELISA Kit	75pg/ml	125-8000pg/ml
EH0158	Human HBEGF(Proheparin-binding EGF-like growth factor) ELISA Kit	18.75pg/ml	31.2-2000pg/ml
EH0157	Human GzmB(Granzyme B) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0156	Human GNLY(Granulysin) ELISA Kit	0.141ng/ml	0.234-15ng/ml
EH0153	Human GM-CSF(Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit	4.688pg/ml	7.813-500pg/ml
EH0152	Human GH(Growth Hormone) ELISA Kit	15pg/ml	25-1600pg/ml
EH0151	Human GDNF(Glial cell line-derived neurotrophic factor) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EH0150	Human GDF15(Growth Differentiation Factor 15) ELISA Kit	14.063pg/ml	23.438-1500pg/ml
EH0148	Human GAL9(Galectin-9) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0145	Human GAL3(Galectin-3) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0138	Human FSTL3(Follistatin-related protein 3) ELISA Kit	0.188ng/ml	0.312-20ng/ml
EH0141	Human Fractalkine(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0134	Human FN1(Fibronectin) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EH0130	Human FGF21(Fibroblast growth factor 21) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EH0302	Human TNF α (Tumor Necrosis Factor Alpha) ELISA Kit	1pg/ml	15.6-1000pg/ml
EH0252	Human PD-1/PDCD1(Programmed Cell Death Protein 1) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0357	Human EPO(Erythropoietin) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EM0004	Mouse CXCL10/IP-10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	18.75pg/ml	31.2-2000pg/ml
EM0006	Mouse CXCL4/PF4(Platelet Factor 4) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EM0010	Mouse MPO(Myeloperoxidase) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EM0020	Mouse BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	2pg/ml	31.2-2000pg/ml
EM0022	Mouse BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EM0050	Mouse CD40L(Cluster of Differentiation 40 Ligand) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EM0057	Mouse Chi3l1(Chitinase-3-like protein 1) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EM0061	Mouse hs-CRP (High Sensitivity C-Reactive Protein) ELISA Kit	0.234ng/ml	0.391-25ng/ml
EM0073	Mouse EPO(Erythropoietin) ELISA Kit	28.125pg/ml	47-3000pg/ml
EM0086	Mouse G-CSF(Granulocyte Colony Stimulating Factor) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0088	Mouse GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EM0089	Mouse GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0093	Mouse IFN γ (Interferon Gamma) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0094	Mouse IGF1(Insulin-like growth factor 1) ELISA Kit	37.5pg/ml	62.4-4000pg/ml
EM0098	Mouse IGFBP-3(Insulin-Like Growth Factor Binding Protein 3) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EM0100	Mouse IL-10(Interleukin-10) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0103	Mouse IL-13(Interleukin 13) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EM0108	Mouse IL-1 α (Interleukin 1 Alpha) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0109	Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit	7.5pg/ml	12.5-800pg/ml
EM0112	Mouse IL-2(Interleukin 2) ELISA Kit	18.75pg/ml	31.25-2000pg/ml

Hot Selling Kits

Cat. No.	Product Name	Sensitivity	Detection Range
EM0113	Mouse IL-21(Interleukin 21) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EM0114	Mouse IL-23(Interleukin 23) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0116	Mouse IL-3(Interleukin 3) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0118	Mouse IL-33(Interleukin 33) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0119	Mouse IL-4(Interleukin 4) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0120	Mouse IL-5(Interleukin 5) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0121	Mouse IL-6(Interleukin 6) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
EM0127	Mouse LN(Laminin) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EM0129	Mouse Leptin ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EM0135	Mouse MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
EM0142	Mouse MMP-2(Matrix Metalloproteinase 2) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EM0148	Mouse NGF/NGF β (Nerve growth factor) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EM0151	Mouse OPN(Osteopontin) ELISA Kit	23.438pg/ml	39.063-2500pg/ml
EM0174	Mouse SDF-1(Stromal Cell Derived Factor 1) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EM0176	Mouse TGF- β 1(Transforming Growth Factor β 1) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EM0178	Mouse TGF- β 3(Transforming Growth Factor β 3) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
EM0183	Mouse TNF α (Tumor Necrosis Factor Alpha) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0008	Rat BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0010	Rat BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
ER0012	Rat IFN- γ (Interferon γ) ELISA Kit	9.4pg/ml	15.6-1000pg/ml
ER0019	Rat ET-1(Endothelin 1) ELISA Kit	4.688pg/ml	7.813-500pg/ml
ER0028	Rat ICAM-1(Intercellular Adhesion Molecule-1) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
ER0029	Rat IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0030	Rat IGF-1(Insulin-Like Growth Factor 1) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
ER0033	Rat IL-10(Interleukin-10) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0034	Rat IL-13(Interleukin-13) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0035	Rat IL-17(Interleukin-17) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0036	Rat IL-18(Interleukin-18) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0039	Rat IL-2(Interleukin-2) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0040	Rat IL-3(Interleukin-3) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0041	Rat IL-4(Interleukin-4) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0042	Rat IL-6(Interleukin-6) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
ER0043	Rat Kim-1(Kidney Injury Molecule 1) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
ER0047	Rat MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0049	Rat MIP-1 α /CCL3(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0051	Rat MMP-2(Matrix Metalloproteinase 2) ELISA Kit	0.188ng/ml	0.313-20ng/ml
ER0054	Rat NGF/NGF β (Beta-nerve growth factor) ELISA Kit	23.438pg/ml	39.063-2500pg/ml
ER0055	Rat NT-3(Neurotrophin 3) ELISA Kit	18.75pg/ml	31.2-2000pg/ml
ER0059	Rat Rantes(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0061	Rat TGF β I(Transforming Growth Factor Beta Induced Protein) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0064	Rat TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	23.438pg/ml	39.063-2500pg/ml
ER0065	Rat TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	0.188ng/ml	0.313-20ng/ml
ER0069	Rat VEGF(Vascular Endothelial cell Growth Factor) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0074	Rat Afp(Alpha-fetoprotein) ELISA Kit	0.188ng/ml	0.312-20ng/ml
ER0080	Rat BMP-4(Bone morphogenetic protein 4) ELISA Kit	9.375pg/ml	15.625-1000pg/ml

Hot Selling Kits

Cat. No.	Product Name	Sensitivity	Detection Range
ER0081	Rat IGFBP-1(Insulin-like growth factor-binding protein 1) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0085	Rat VEGF-B(Vascular Endothelial Cell Growth Factor B) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0087	Rat Flt1/VEGFR1(Vascular endothelial growth factor receptor 1) ELISA Kit	0.094ng/ml	0.156-10ng/ml
ER0093	Rat IL-11(Interleukin-11) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0094	Rat Ace(Angiotensin-converting enzyme) ELISA Kit	0.188ng/ml	0.313-20ng/ml
ER0107	Rat IL-15(Interleukin-15) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0114	Rat IL-5(Interleukin-5) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
ER0115	Rat Lep(Leptin) ELISA Kit	0.094ng/ml	0.156-10ng/ml
ER0124	Rat SELP(P-selectin) ELISA Kit	0.188ng/ml	0.313-20ng/ml
ER0134	Rat Pcsk9(Proprotein convertase subtilisin/kexin type 9) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0139	Rat MMP-9(Matrix Metalloproteinase-9) ELISA Kit	0.469ng/ml	0.78-50ng/ml
ER0142	Rat MPO(Myeloperoxidase) ELISA Kit	0.469ng/ml	0.78-50ng/ml
ER0143	Rat Casp3(Caspase 3) ELISA Kit	0.188ng/ml	0.313-20ng/ml
ER0145	Rat PDGFB(Platelet-derived growth factor subunit B) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ER0146	Rat IL-7(Interleukin-7) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
EMK0019	Monkey Cys-C(Cystatin C) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EMK0021	Monkey EPO(Erythropoietin) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EMK0028	Monkey IFN- α (Interferon Alpha) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EMK0029	Monkey IFN- β (Interferon Beta) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EMK0030	Monkey IFN- γ (Interferon Gamma) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EMK0033	Monkey IGF-1(Insulin Like Growth Factor 1) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EMK0034	Monkey IgG(Immunoglobulin G) ELISA Kit	0.938ug/ml	1.563-100ug/ml
EMK0036	Monkey IL-10(Interleukin-10) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EMK0039	Monkey IL-2(Interleukin 2) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EU0253	ST/5-HT(5-hydroxytryptamine) ELISA Kit	9.375ng/ml	15.625-1000ng/ml
EU0254	Androstenedione(Androstenedione) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EU0380	Pregnenolone(Pregnenolone) ELISA Kit	9.375pg/ml	15.6-1000pg/ml
EU0381	FA/VB9(Folic Acid/Vitamin B9) ELISA Kit	46.9pg/ml	78-5000pg/ml
EU0390	E2(Estradiol) ELISA Kit	0.422ng/ml	0.703-45ng/ml
EU0391	COR(Cortisol) ELISA Kit	0.234ng/ml	0.391-25ng/ml
EU0392	DA(Dopamine) ELISA Kit	0.938ng/ml	1.56-100ng/ml
EU0398	PG(Progesterone) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EU0400	T(Testosterone) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EG0036	Goat TGF β 1 ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EGP0051	Guinea Pig TGF- β 1(Transforming Growth Factor Beta 1) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EB0010	Bovine TGF- β (Transforming Growth Factor Beta) ELISA Kit	18.75pg/ml	31.2-2000pg/ml
ECM0013	Camel IgG(Immunoglobulin G) ELISA Kit	1.875ng/ml	3.13-200ng/ml
ECA0019	Canine TGF- β 1(Transforming Growth Factor β 1) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EGP0020	Guinea pig IgG(Immunoglobulin G) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EHA0009	Hamster TGF- β (Transforming Growth Factor beta) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EP0154	Porcine TGF- β 1(Transforming Growth Factor β 1) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
ESH0060	Sheep IgG(Immunoglobulin G) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EB0148	Bovine serum albumin (BSA) ELISA Kit	0.094ug/ml	0.156-10ug/ml
ECA0011	Canine IgG(Immunoglobulin G) ELISA Kit	0.469ng/ml	0.781-50ng/ml
ERB0119	Rabbit TGF- β 1(Transforming Growth Factor β 1) ELISA Kit	18.75pg/ml	31.25-2000pg/ml

Hot Selling Kits

Cat. No.	Product Name	Sensitivity	Detection Range
EH0029	Human ADAM12(Disintegrin and metalloproteinase domain-containing protein 12) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0359	Human AFP(alpha Fetoprotein) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH1268	Human AFP-L3(Alpha-fetoprotein Lens culinaris agglutinin 3) ELISA Kit	0.375ng/ml	0.625-40ng/ml
EM1856	Mouse AFP-L3(Alpha-Fetoprotein L3) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH0424	Human B2M(Beta2-microglobulin) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EM0209	Mouse B2m(Beta-2-microglobulin) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EH0669	Human BAX(Apoptosis regulator BAX) ELISA Kit	0.469ng/ml	0.781-50ng/ml
EH0658	Human Bcl-2(B-cell Leukemia/Lymphoma 2) ELISA Kit	0.938ng/ml	1.56-100ng/ml
EH1088	Human BRCA1(Breast cancer type 1 susceptibility protein) ELISA Kit	0.188ng/ml	0.312-20ng/ml
EH1532	Human BRCA2(Breast cancer type 2 susceptibility protein) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH0361	Human CA-125(Carbohydrate Antigen 125) ELISA Kit	0.938U/ml	1.563-100U/ml
EM0887	Mouse CA125(Carbohydrate Antigen 125) ELISA Kit	1.875U/ml	3.125-200U/ml
ER0788	Rat CA125(Carbohydrate Antigen 125) ELISA Kit	4.688U/ml	7.813-500U/ml
EH0363	Human CA19-9(Carbohydrate antigen19-9) ELISA Kit	1.875U/ml	3.125-200U/ml
EH0595	Human CASP1(Caspase-1) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EH0073	Human CD147(Extracellular Matrix Metalloproteinase Inducer) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0634	Human CD276(CD276 antigen) ELISA Kit	0.469ng/ml	0.781-50ng/ml
EH0090	Human CEA(Carcinoembryonic Antigen) ELISA Kit	0.188ng/ml	0.312-20ng/ml
EM0916	Mouse CEA(Carcinoembryonic Antigen) ELISA Kit	37.5pg/ml	62.5-4000pg/ml
EH0486	Human CGA(Chromogranin A) ELISA Kit	0.469ng/ml	0.781-50ng/ml
EM0924	Mouse CHGA(Chromogranin A) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EH1342	Human CLDN3(Claudin-3) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EM0538	Mouse Cldn3(Claudin-3) ELISA Kit	0.094ng/ml	0.156-10ng/ml
ER0457	Rat Cldn3(Claudin-3) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0094	Human c-MET(Hepatocyte growth factor receptor) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EH1356	Human CNDP1(Beta-Ala-His dipeptidase) ELISA Kit	0.938ng/ml	1.56-100ng/ml
ER1663	Rat CNDP1(Beta-Ala-His dipeptidase) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EH0759	Human CRYAB(Alpha-crystallin B chain) ELISA Kit	0.094ng/ml	0.156-10ng/ml
EH0377	Human CT(Calcitonin) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EM0954	Mouse CT(Calcitonin) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
ER0867	Rat CT(Calcitonin) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH0364	Human CYFRA21-1(Cytokeratin Fragment Antigen 21-1) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH2114	Human EDN2(Endothelin-2) ELISA Kit	0.469pg/ml	0.781-50pg/ml
EM0727	Mouse Edn2(Endothelin-2) ELISA Kit	46.9pg/ml	78-5000pg/ml
EH0010	Human EGFR(Epidermal growth factor receptor) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EM1602	Mouse EGFR(Epidermal Growth Factor Receptor) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
ER0915	Rat EGFR(Epidermal Growth Factor Receptor) ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EH2381	Human EPX(Eosinophil peroxidase) ELISA Kit	0.938ng/ml	1.56-100ng/ml
EH4301	Human Total ER(Estrogen Receptor)ELISA Kit	0.281ng/ml	0.469-30ng/ml
EH2013	Human EREG(Proepiregulin) ELISA Kit	9.375pg/ml	15.625-1000pg/ml
EH3014	Human ER α (Estrogen Receptor Alpha) ELISA Kit	0.375ng/ml	0.625-40ng/ml
EM1007	Mouse ER α (Estrogen Receptor Alpha) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EH3015	Human ER β (Estrogen Receptor Beta) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EM1008	Mouse ER β (Estrogen Receptor Beta) ELISA Kit	0.188ng/ml	0.313-20ng/ml
ER0424	Rat Esr1(Estrogen receptor) ELISA Kit	0.188ng/ml	0.312-20ng/ml

Citation Selection



Cat.No.	Product name	Journal	IF
EH4233	Human PGE2(Prostaglandin E2) ELISA Kit	Cell	31.398
EH9300	Human IFITM2(Interferon-induced transmembrane protein 2)ELISA Kit	Hepatology	14.079
EH0302	Human TNF α (Tumor Necrosis Factor Alpha) ELISA Kit	Journal of Crohn's and Colitis	6.992
EH0287	Human TGF β 1(Transforming Growth Factor Beta 1) ELISA Kit	Journal of Crohn's and Colitis	6.992
EM0144	Mouse MMP-9(Matrix Metalloproteinase 9) ELISA Kit	Journal of Experimental & Clinical Cancer Research	6.217
EM0363	Mouse Ctsk(Cathepsin K) ELISA Kit	Journal of Experimental & Clinical Cancer Research	6.217
ER1901	Rat GLUR2(Glutamate receptor 2) Elisa Kit	Molecular Neurobiology	6.19
EM0310	Mouse HIF-1 α (Hypoxia-inducible factor 1-alpha) ELISA Kit	The FASEB Journal	5.5
ER1393	Rat TNF- α (Tumor Necrosis Factor Alpha)ELISA Kit	Cellular Physiology and Biochemistry	5.5
EH2862	Human COL15(Collagen Type XV) ELISA Kit	Cellular Physiology and Biochemistry	5.5
ER1619	Rat OXTR(Oxytocin receptor) ELISA Kit	FASEB Journal	5.498
EU3120	6-hydroxymelatonin sulfate ELISA Kit	Clinical Nutrition	5.496
ER0094	Rat Ace(Angiotensin-converting enzyme) ELISA Kit	British Journal of Pharmacology	5.49
ER0609	Rat Ace2(Angiotensin-converting enzyme 2) ELISA Kit	British Journal of Pharmacology	5.49
ER1637	Rat ANG II (Angiotensin II) ELISA Kit	British Journal of Pharmacology	5.49
ER1616	Rat Ang1-7(Angiotensin 1-7) ELISA Kit	British Journal of Pharmacology	5.49
ECH0069	Chicken cTnI(Cardiac Troponin-I) ELISA Kit	Environmental Pollution	5.291
EH3732	Human RTN4(Reticulon 4) ELISA Kit	Oncotarget	5.168

Citation Selection

Cat.No.	Product name	Journal	IF
EH0009	Human EGF(Epidermal Growth Factor) ELISA Kit	Oncotarget	5.168
EH0327	Human VEGF(Vascular Endothelial Cell Growth Factor) ELISA Kit	Oncotarget	5.168
EM0422	Mouse ApoA1(Apolipoprotein A-I) ELISA Kit	Molecular Nutrition & Food Research	5.151
ER0332	Rat SOD1(Superoxide dismutase [Cu-Zn])ELISA Kit	Molecular Nutrition & Food Research	5.151
ER0870	Rat Tn-I(Troponin I) ELISA Kit	Clinical Science	4.996
EU2543	VD3(Vitamin D3) ELISA Kit	Oxidative Medicine and Cellular Longevity	4.936
EH3057	Human FG(Fibrinogen) ELISA Kit	Cerebral Blood Flow & Metabolism	4.929
ER0049	Rat MIP-1 α /CCL3 ELISA Kit	Scientific Reports	4.847
EH3057	Human FG(Fibrinogen) ELISA Kit	Scientific Reports	4.847
EU2571	LXA4(Lipoxin A4) ELISA Kit	Journal of Dental Research	4.755
EM1503	Mouse PGE2(Prostaglandin E2) ELISA Kit	Journal of Dental Research	4.755
EU0198	AP(Allopregnanolone) ELISA Kit	The Journal of Pain	4.519
ER1384	Rat TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	Oxidative Medicine and Cellular Longevity	4.492
EM0109	Mouse IL-1b(Interleukin 1 Beta) ELISA Kit	Molecular and Cellular Endocrinology	4.4
EM0260	Mouse INS(Insulin) ELISA Kit	Scientific Reports	4.26
ER1138	Rat MAPK (P38 Mitogen-Activated Protein Kinase) ELISA Kit	Scientific Reports	4.26
EM0572	Mouse O3far1(Omega-3 fatty acid receptor 1) ELISA Kit	Arthritis Research & Therapy	4.2
EM0142	Mouse MMP-2(Matrix Metalloproteinase 2) ELISA Kit	Scientific Reports	4.122
EU0253	5-HT(5-hydroxytryptamine) ELISA Kit	Scientific Reports	4.122
EM0010	Mouse MPO(Myeloperoxidase) ELISA Kit	Scientific Reports	4.122
EM0864	Mouse A β 42(Amyloid Beta 42) ELISA Kit	Toxicological Sciences	4.081
EM0121	Mouse IL-6(Interleukin 6) ELISA Kit	Cellular Physiology	4.08
ER1113	Rat INS(Insulin) ELISA Kit	Frontiers in Microbiology	4.019